

Use

To use pHIS2 in a one-hybrid assay, clone one or more copies of a *cis*-acting DNA target sequence into the MCS. Then introduce the plasmid into competent yeast cells using the transformation protocol in the Matchmaker Library Construction & Screening Kits User Manual (PT3529-1). In contrast to the original Matchmaker One-Hybrid System, this reporter vector does not need to be integrated into the yeast genome. Instead, it is maintained as an episome throughout the assay. Inserting your target element may alter the level of background *HIS3* expression. Therefore, constructs should be tested for background (leaky) *HIS3* expression before you start a one-hybrid analysis. Background growth due to leaky *HIS3* expression is controlled by adding 3-AT to the selection medium, as described in the User Manual (PT3529-1).

Location of Features

- Multiple cloning sites: 1–41
- *HIS3* gene: 152–811
- Fragment containing the *HIS3* 3' UTR & Termination sequence: 812–1446
- *TRP1* gene: 2855–3529
- Fragment containing the Col E1 *E. coli* origin of replication: 4165–4615
- Kanamycin-resistance gene: 5605–4811
- *CEN6/ARS4* sequences: 6254–5737

Propagation in *E. coli*

- Suitable host strains: DH5 α and other general purpose strains.
- Selectable marker: plasmid confers resistance to kanamycin (50 μ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: Col E1
- Copy number: low

Propagation in *S. cerevisiae*

- Suitable host strains: AH109(*MATa*), Y187(*MAT α*), Y190(*MATa*), SFY526(*MATa*), CG1945(*MATa*), HF7c(*MATa*)
- Selectable marker: *TRP1*
- *S. cerevisiae* origin: *CEN6/ARS4*

References

1. Sikorski, R. S. & Hieter, P. (1989) *Genetics* **122**:19–27.
2. Rose, M. D. & Broach, J. R. (1991) *Methods Enzymol.* **194**:195–230.

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech Laboratories, Inc.. This vector has not been completely sequenced.

Notice to Purchaser

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes, nor is it intended for human use. Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Clontech Laboratories, Inc.

Clontech, Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.
Clontech is a Takara Bio Company. ©2005