

pcDNA[™] 3.1/*myc*-His(-) A, B, and C

Catalog no. V855-20

Version G
October 28, 2010
28-0128



www.invitrogen.com
tech_service@invitrogen.com

Table of Contents

Table of Contents.....	iii
Important Information.....	iv
Overview.....	1
pcDNA™ 3.1/ <i>myc</i> -His(-) A, B, and C Vectors.....	2
Cloning into pcDNA™ 3.1/ <i>myc</i> -His(-) A, B, and C.....	4
General Guidelines and Special Information.....	8
pcDNA™ 3.1/ <i>myc</i> -His(-)/ <i>lacZ</i>	10
Technical Service.....	11
Purchaser Notification.....	13
References.....	14

Important Information

Contents 20 µg each pcDNA[™]3.1/*myc*-His(-) A, B, and C, lyophilized
20 µg pcDNA[™]3.1/*myc*-His(-)/*lacZ*, lyophilized

Shipping/Storage Lyophilized plasmids are shipped at room temperature and should be stored at -20°C.

Product Qualification The pcDNA[™]3.1/*myc*-His(-) and pcDNA[™]3.1/*myc*-His(-)/*lacZ* vectors are qualified by restriction enzyme digestion with specific restriction enzymes as listed below. Restriction digests must demonstrate the correct banding pattern when electrophoresed on an agarose gel (see below).

Vector	Restriction Enzymes	Expected Results (bp)
pcDNA [™] 3.1/ <i>myc</i> -His(-) A	<i>Apa</i> I	92, 5430
	<i>Sna</i> B I	5522
	<i>Xba</i> I	5522
pcDNA [™] 3.1/ <i>myc</i> -His(-) B	<i>Apa</i> I	5520
	<i>Sna</i> B I	5520
	<i>Xba</i> I	86, 5434
pcDNA [™] 3.1/ <i>myc</i> -His(-) C	<i>Apa</i> I	5521
	<i>Sna</i> B I	412, 5109
	<i>Xba</i> I	5521
pcDNA [™] 3.1/ <i>myc</i> -His(-)/ <i>lacZ</i>	<i>Nhe</i> I	8592
	<i>Xba</i> I	8592
	<i>Xba</i> I + <i>Eco</i> R V	1151, 7441

Overview

Introduction

pcDNA[™] 3.1/*myc*-His(-) A, B, and C are 5.5 kb vectors designed for overproduction of recombinant proteins in mammalian cell lines. The vectors are supplied in three reading frames to facilitate in frame cloning, with a C-terminal peptide containing a polyhistidine metal-binding tag and the *myc* (*c-myc*) epitope. The human cytomegalovirus immediate-early (CMV) promoter provides high-level expression in a wide range of mammalian cells. In addition, the vector will replicate episomally in cell lines that are latently infected with SV40 or that express the SV40 large T antigen (e.g. COS7). High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The control plasmid, pcDNA[™] 3.1/*myc*-His(-)/*lacZ*, is the pcDNA[™] 3.1/*myc*-His(-) A vector with a 3.2 kb fragment containing the β -galactosidase gene cloned in frame with the C-terminal peptide (see page 10). It is included for use as a positive control for transfection, expression, and detection in the cell line of choice.

Experimental Outline

Use the following outline to clone and express your gene of interest in pcDNA[™] 3.1/*myc*-His(-).

- Consult the multiple cloning sites described on pages 5-7 to determine which vector (A, B, or C) should be used to clone your gene in frame with the C-terminal *myc* epitope and the polyhistidine tag.
- Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on 50 to 100 μ g/ml ampicillin.
- Analyze your transformants for the presence of insert by restriction digestion.
- Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in frame with the C-terminal peptide.
- Transfect your construct into the cell line of choice using your own method of transfection. Generate a stable cell line, if desired.
- Test for expression of your recombinant gene by western blot analysis or functional assay. If you do not have an antibody to your protein, you may use the Anti-*myc* Antibody (Catalog no. R950-25) or the Anti-Hi (C-term) Antibody (Catalog no. R930-25) to detect your recombinant protein.
- To purify your recombinant protein, you may use metal-chelating resin such as ProBond[™]. ProBond[™] resin is available separately (see below for ordering information).

ProBond[™] Resin

Ordering information for ProBond[™] resin is provided below.

Item	Amount	Catalog no.
ProBond [™] Purification System	6 x 2 ml precharged ProBond [™] resin columns and buffers for native and denaturing purification	K850-01
ProBond [™] Purification System with Anti- <i>myc</i> -HRP Antibody	1 Kit	K852-01
ProBond [™]	50 ml	R801-01
	150 ml	R801-15

pcDNA™ 3.1/myc-His(-) A, B, and C Vectors

Features of pcDNA™ 3.1/myc- His(-)

pcDNA™ 3.1/myc-His(-) A (5522 bp), pcDNA™ 3.1/myc-His(-) B (5520 bp), and pcDNA™ 3.1/myc-His(-) C (5521 bp) contain the following elements. All features have been functionally tested.


Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Allows efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
Multiple cloning site in three reading frames	Allows insertion of your gene and facilitates cloning in frame with the C-terminal polyhistidine tag
<i>myc</i> epitope (<i>c-myc</i>)	Allows detection of your recombinant protein with the Anti- <i>myc</i> Antibody (Catalog no. R950-25) or Anti- <i>myc</i> -HRP Antibody (Catalog no. R951-25) (Evan <i>et al.</i> , 1985)
C-terminal polyhistidine tag	Allows purification of your recombinant protein on metal-chelating resin such as ProBond™ In addition, the C-terminal polyhistidine tag is the epitope for the Anti-His (C-term) Antibody (Catalog no. R930-25) and the Anti-His (C-term)-HRP Antibody (Catalog no. R931-25)
BGH reverse priming site	Allows sequencing through the insert
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992)
f1 origin	Allows rescue of single-stranded DNA
SV40 early promoter and origin	Allows efficient, high-level expression of the neomycin resistance gene and episomal replication in cells expressing SV40 large T antigen
Neomycin (Geneticin®) resistance gene	Selection of stable transfectants in mammalian cells (Southern and Berg, 1982)
SV40 polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA
pUC origin	High-copy number replication and growth in <i>E. coli</i>
Ampicillin resistance gene (β -lactamase)	Selection of vector in <i>E. coli</i>

continued on next page

pcDNA™ 3.1/*myc*-His(-) A, B, and C Vectors, continued

Map of pcDNA™ 3.1/*myc*-His(-)

The figure below summarizes the features of the pcDNA™ 3.1/*myc*-His(-) vectors. The sequences for pcDNA™ 3.1/*myc*-His(-) A, B, and C are available for downloading from our Web site (www.invitrogen.com) or from Technical Service (see page 11). Details of the multiple cloning sites for pcDNA™ 3.1/*myc*-His(-) A, B, and C are shown on pages 5-7.



**Comments for pcDNA™ 3.1/*myc*-His(-) A:
5522 nucleotides**

- CMV promoter: bases 209-863
- T7 promoter/priming site: bases 863-882
- Multiple cloning site: bases 895-1006
- myc* epitope: bases 1007-1036
- Polyhistidine tag: bases 1052-1069
- BGH reverse priming site: bases 1113-1130
- BGH polyadenylation signal: bases 1116-1343
- f1 origin: bases 1389-1817
- SV40 promoter and origin: bases 1844-2152
- Neomycin resistance gene: bases 2227-3021
- SV40 polyadenylation signal: bases 3195-3325
- pUC origin: bases 3708-4381
- Ampicillin resistance gene: bases 4526-5386 (complementary strand)

* There are two *Apa* I sites in version A only.

** There are two *Xba* I sites in version B only.

Cloning into pcDNA™ 3.1/myc-His(-) A, B, and C

General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

E. coli Strain

Many *E. coli* strains are suitable for the growth of this vector including TOP10F'. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*) and endonuclease A-deficient (*endA*).

For your convenience, TOP10F' is available as chemically competent or electrocompetent cells from Invitrogen.

Item	Quantity	Catalog no.
Electrocomp™ TOP10F'	5 x 80 µl	C665-55
One Shot® TOP10F' (chemically competent cells)	21 x 50 µl	C3030-03

Maintenance of pcDNA™ 3.1/myc-His(-)

To propagate and maintain the pcDNA™ 3.1/myc-His(-) vectors, we recommend resuspending each vector in 20 µl sterile water to make a 1 µg/µl stock solution. Store the stock solution at -20°C.

Use this stock solution to transform a *recA*, *endA* *E. coli* strain like TOP10, TOP10F', DH5α, JM109, or equivalent. Select transformants on LB plates containing 50-100 µg/ml ampicillin.

Cloning Considerations

Your insert should contain a Kozak translation initiation sequence for proper initiation of translation (Kozak, 1987; Kozak, 1990). An example of a Kozak consensus sequence is provided below. Note that other sequences are possible, but the G or A at position -3 and the G at position +4 are the most critical nucleotides for function (shown in bold). The ATG initiation codon is shown underlined.

(G/A)NN**ATGG**

To express your gene as a recombinant fusion protein, you must clone your gene in frame with the C-terminal peptide. The vector is supplied in three reading frames to facilitate cloning. See pages 5-7 to develop a cloning strategy.

If you wish to express your protein **without** the C-terminal peptide, be sure to include a stop codon.

continued on next page

Cloning into pcDNA™ 3.1/myc-His(-) A, B, and C, continued

Multiple Cloning Site of Version A

Below is the multiple cloning site for pcDNA™ 3.1/myc-His(-) A. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. **Note that there is a stop codon after the Xba I site.** The multiple cloning site has been confirmed by sequencing and functional testing. **The complete sequence of pcDNA™ 3.1/myc-His(-) A is available for downloading from our Web site (www.invitrogen.com) or from Technical Service (see page 11).**

```

741 AAATGTCGTA ACAACTCCGC CCCATTGACG CAAT CAAATGGGCG GTAGGCGTGT ACGGTGGGAG
      TATA                                     putative transcriptional start
801 GTCTATATAA GCAGAGCTCT CTGGCTAACT AGAGAACCCA CTGCTTACTG GCTTATCGAA
      T7 promoter/priming site                Nhe I                               Apa I* Xba I
861 ATTAATACGA CTCACTATAG GGAGACCCAA GCTGG CTA GCG TTT AAA CGG GCC CTC
      Leu Ala Phe Lys Arg Ala Leu
      Xho I   Not I           BstX I*   EcoR V       EcoR I           BstX I*
917 TAG ACT CGA GCG GCC GCC ACT GTG CTG GAT ATC TGC AGA ATT CCA CCA CAC
   *** Thr Arg Ala Ala Ala Thr Val Leu Asp Ile Cys Arg Ile Pro Pro His
      BamH I           Asp718 I   Kpn I   Hind III       Apa I*       myc epitope
968 TGG ACT AGT GGA TCC GAG CTC GGT ACC AAG CTT GGG CCC GAA CAA AAA CTC
   Trp Thr Ser Gly Ser Glu Leu Gly Thr Lys Leu Gly Pro Glu Gln Lys Leu
      Polyhistidine tag
1019 ATC TCA GAA GAG GAT CTG AAT AGC GCC GTC GAC CAT CAT CAT CAT CAT CAT
   Ile Ser Glu Glu Asp Leu Asn Ser Ala Val Asp His His His His His His
      Afl II
1070 TGA GTTTAAACGG TCTCCAGCTT AAGTTTAAAC CGCTGATCAG CCTCGACTGT GCCTTCTA
   ***
  
```

*Please note that there are two *BstX I* sites and two *Apa I* sites in the polylinker.

continued on next page

Cloning into pcDNA™ 3.1/myc-His(-) A, B, and C, continued

Multiple Cloning Site of Version B

Below is the multiple cloning site for pcDNA™ 3.1/myc-His(-) B. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. **Note that there is a stop codon after the *Nhe* I site.** The multiple cloning site has been confirmed by sequencing and functional testing. **The complete sequence of pcDNA™ 3.1/myc-His(-) B is available for downloading from our Web site (www.invitrogen.com) or from Technical Service (see page 11).**

```

741 AAATGTCGTA ACAACTCCGC CCCATTGACG CAAATGGGCG GTAGGCGTGT ACGGTGGGAG
      TATA
801 GTCTATATAA GCAGAGCTCT CTGGCTAACT AGAGAACCCA CTGCTTACTG GCTTATCGAA
      T7 promoter/priming site
861 ATTAATACGA CTCACTATAG GGAGACCCAA GCT GGC TAG CGT TTA AAC GGG CCC
      Gly *** Arg Leu Asn Gly Pro
      Nhe I
      Apa I
915 TCT AGA CTC GAG CGG CCG CCA CTG TGC TGG ATA TCT GCA GAA TTC CAC
      Ser Arg Leu Glu Arg Pro Pro Leu Cys Trp Ile Ser Ala Glu Phe His
      Xba I* Xho I Not I BstX I* EcoR V EcoR I
963 CAC ACT GGA CTA GTG GAT CCG AGC TCG GTA CCA AGC TTT CTA GAA CAA
      His Thr Gly Leu Val Asp Pro Ser Ser Val Pro Ser Phe Leu Glu Gln
      BstX I* BamH I Asp718 I Kpn I Hind III Xba I* myc epitope
1011 AAA CTC ATC TCA GAA GAG GAT CTG AAT AGC GCC GTC GAC CAT CAT CAT
      Lys Leu Ile Ser Glu Glu Asp Leu Asn Ser Ala Val Asp His His His
      Polyhistidine tag
1059 CAT CAT CAT TGA GTTTAAACGG TCTCCAGCTT AAGTTTAAAC CGCTGATCAG
      His His His ***
      Afl II
1111 CCTCGACTGT GCCTTCTAGT TGCCAGCCAT CTGTTGTTTG CCCCTCCCC GTGCCTTCCT
      BGH Reverse priming site
  
```

*Please note that there are two *Bst*X I sites and two *Xba* I sites in the polylinker.

continued on next page

Cloning into pcDNA™ 3.1/myc-His(-) A, B, and C, continued

Multiple Cloning Site of Version C

Below is the multiple cloning site for pcDNA™ 3.1/myc-His(-) C. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. **Note that there are two stop codons. One is located between the *Nhe* I site and the *Apa* I site, and the other is located before the *Bam*H I site.** This means the 3' cloning site must be either *Bam*H I, *Kpn* I (*Asp*718 I), or *Hind* III if you wish to clone your gene into pcDNA™ 3.1/myc-His(-) C so that it is expressed without interruption. The multiple cloning site has been confirmed by sequencing and functional testing. **The complete sequence of pcDNA™ 3.1/myc-His(-) C is available for downloading from our Web site (www.invitrogen.com) or from Technical Service (see page 11).**

```

741 AAATGTCGTA ACAACTCCGC CCCATTGACG CAAT
      |
      | CAAT
      |
      | CAAATGGGCG GTAGGCGTGT ACGGTGGGAG
      |
      |
801 TATA
      |
      | putative transcriptional start
      |
      | GTCTATATAA GCAGAGCTCT CTGGCTAACT AGAGAACCCA CTGCTTACTG GCTTATCGAA
      |
      |
861 T7 promoter/priming site
      |
      | ATTAATACGA CTCACTATAG GGAGACCCAA GCTG Nhe I
      |
      | GCT AGC GTT TAA ACG GGC CCT Apa I Xba I
      |
      | Ala Ser Val *** Thr Gly Pro
      |
      |
916 Xho I Not I BstX I* EcoR V EcoR I
      |
      | CTA GAC TCG AGC GGC CGC CAC TGT GCT GGA TAT CTG CAG AAT TCC ACC
      |
      | Leu Asp Ser Ser Gly Arg His Cys Ala Gly Tyr Leu Gln Asn Ser Thr
      |
      |
964 BstX I* BamH I Asp718 I Kpn I Hind III
      |
      | ACA CTG GAC TAG TGG ATC CGA GCT CGG TAC CAA GCT TAC GTA GAA CAA
      |
      | Thr Leu Asp *** Trp Ile Arg Ala Arg Tyr Gln Ala Tyr Val Glu Gln
      |
      |
1012 myc epitope
      |
      | AAA CTC ATC TCA GAA GAG GAT CTG AAT AGC GCC GTC GAC CAT CAT CAT
      |
      | Lys Leu Ile Ser Glu Glu Asp Leu Asn Ser Ala Val Asp His His His
      |
      |
1060 Afl II
      |
      | CAT CAT CAT TGA GTTTAAACG GTCTCCAGCT TAAGTTTAAA CCGCTGATCA
      |
      | His His His ***
      |
      |
1111 BGH Reverse priming site
      |
      | GCCTCGACTG TGCCTTCTAG TTGCCAGCCA TCTGTTGTTT GCCCCTCCCC CGTGCCTTCC
  
```

*Please note that there are two *Bst*X I sites in the polylinker.

General Guidelines and Special Information

Introduction

The following guidelines and recommendations are provided for your convenience. If you need more details about the techniques discussed, refer to the general molecular biology references in the **Reference** section.

E. coli Transformation

Transform your ligation mixtures into a competent *recA*, *endA* *E. coli* strain (e.g. TOP10, TOP10F', DH5 α) and select on LB plates containing 50-100 μ g/ml ampicillin. Select 10-20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct to confirm that your gene is fused in frame with the *myc* epitope and the C-terminal polyhistidine tag. We suggest using the T7 Promoter and BGH Reverse primer sequences. Refer to the diagrams on pages 5-7 for the sequence and location of the primer binding sites.

For your convenience, Invitrogen offers a custom primer synthesis service. For more information, refer to our Web site (www.invitrogen.com) or contact Technical Service (see page 11).

Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells and salt will interfere with lipid complexing decreasing transfection efficiency. We recommend isolating DNA using the S.N.A.P.[™] Miniprep Kit (10-15 μ g, Catalog no. K1900-01), the S.N.A.P.[™] MidiPrep Kit (10-200 μ g, Catalog no. K1910-01), or CsCl gradient centrifugation.

Methods of Transfection

For established cell lines (e.g. HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. It is recommended that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology*.

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner and Ringold, 1989) and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). Invitrogen offers the Calcium Phosphate Transfection Kit and a large selection of reagents for transfection. For more information on the reagents available, visit our Web site (www.invitrogen.com) or call Technical Service (see page 11).

Positive Control

pcDNA[™] 3.1/*myc*-His(-)/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 10) and may be used to optimize transfection conditions for your cell line. The gene encoding β -galactosidase is expressed in mammalian cells under the CMV promoter. A successful transfection will result in β -galactosidase expression that can be easily assayed (see next page).

continued on next page

General Guidelines and Special Information, continued

Assay for β -galactosidase Activity

You may assay for β -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the β -Gal Assay Kit (Catalog no. K1455-01) and the β -Gal Staining Kit (Catalog no. K1465-01) for fast, easy detection of β -galactosidase expression.

Geneticin[®] Selective Antibiotic

Geneticin[®] blocks protein synthesis in mammalian cells by interfering with ribosomal function. It is an aminoglycoside, similar in structure to neomycin, gentamycin, and kanamycin. Expression of the bacterial aminoglycoside phosphotransferase gene (APH), derived from Tn5, in mammalian cells results in detoxification of Geneticin[®] (Southern and Berg, 1982).

Geneticin[®] Selection Guidelines

Geneticin[®] is available from Invitrogen (Catalog no. 10131-035). Use as follows:

- Prepare Geneticin[®] in a buffered solution (e.g. 100 mM HEPES, pH 7.3).
- Use 100 to 1000 μ g/ml of Geneticin[®] in complete medium.
- Calculate concentration based on the amount of active drug (check the lot label).
- Test varying concentrations of Geneticin[®] on your cell line to determine the concentration that kills your cells (kill curve). Cells differ in their susceptibility to Geneticin[®].

Cells will divide once or twice in the presence of lethal doses of Geneticin[®], so the effects of the drug take several days to become apparent. Complete selection can take from 3 to 6 weeks of growth in selective medium.

Preparation of Cells for Purification

Use the procedure below to purify recombinant protein from a stable cell line. You will need 5×10^6 to 1×10^7 cells for purification on a 2 ml ProBond[™] column (see ProBond[™] Purification System manual).

1. Seed cells in either five T-75 flasks or 2 to 3 T-175 flasks.
 2. Grow the cells in selective medium until they are 80-90% confluent.
 3. Harvest the cells by treating with trypsin-EDTA for 2 to 5 minutes or by scraping the cells in PBS.
 4. Inactivate the trypsin by diluting in fresh medium (if necessary), and transfer the cells to a sterile microcentrifuge tube.
 5. Centrifuge the cells at 1500 rpm for 5 minutes. You may lyse the cells immediately or freeze in liquid nitrogen and store at -70°C until needed.
-

Lysis of Cells

If you are using ProBond[™] resin, refer to the ProBond[™] Purification System manual for details about sample preparation for chromatography.

If you are using other resin, refer to the manufacturer's instruction for recommendations on sample preparation.



Note

The C-terminal peptide containing the *myc* epitope and the polyhistidine tag will add approximately 3 kDa to the size of your protein.

The size of the LacZ/*myc*-His fusion protein is approximately 121 kDa.


pcDNA™ 3.1/myc-His(-)/lacZ

Description

pcDNA™ 3.1/myc-His(-)/lacZ is a 8592 bp control vector containing the gene for β-galactosidase. pcDNA™ 3.1/myc-His(-) A was digested with *EcoR* V. A 3.2 kb blunt *Sfu* I-*Not* I fragment containing the β-galactosidase gene was then ligated into pcDNA™ 3.1/myc-His(-) A in frame with the C-terminal peptide.

Map of Control Vector

The figure below summarizes the features of the pcDNA™ 3.1/myc-His(-)/lacZ vector. The complete nucleotide sequence for pcDNA™ 3.1/myc-His(-)/lacZ is available for downloading from our Web site (www.invitrogen.com) or from Technical Service (see next page).



Comments for pcDNA™ 3.1/myc-His(-)/lacZ: 8592 nucleotides

CMV promoter: bases 209-863

T7 promoter/priming site: bases 863-882

LacZ ORF: bases 954-4010

myc epitope: bases 4077-4106

Polyhistidine tag: bases 4122-4139

BGH reverse priming site: bases 4183-4200

BGH polyadenylation signal: bases 4186-4413

f1 origin: bases 4459-4887

SV40 promoter and origin: bases 4914-5222

Neomycin resistance gene: bases 5297-6091

SV40 polyadenylation signal: bases 6265-6395

pUC origin: bases 6778-7451

Ampicillin resistance gene: bases 75967-8456 (complementary strand)

* The two *Apa* I sites flanking the *lacZ* gene are unique.

Technical Service

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

www.invitrogen.com

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail:
tech_service@invitrogen.com

Japanese Headquarters:

Invitrogen Japan K.K.
Nihonbashi Hama-Cho Park Bldg. 4F
2-35-4, Hama-Cho, Nihonbashi
Tel: 81 3 3663 7972
Fax: 81 3 3663 8242
E-mail: jpinfo@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

MSDS Requests

To request an MSDS, visit our Web site at www.invitrogen.com. On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

continued on next page

Technical Service, continued

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

Introduction

Use of the pcDNA™3.1/*myc*-His(-) vectors is covered under a number of different licenses including those detailed below.

Limited Use Label License No. 22: Vectors and Clones Encoding Histidine Hexamer

This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nucleic Acids Res.* *15*, 1311-1326.
- Evan, G. I., Lewis, G. K., Ramsay, G., and Bishop, V. M. (1985). Isolation of Monoclonal Antibodies Specific for *c-myc* Proto-oncogene Product. *Mol. Cell. Biol.* *5*, 3610-3616.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Felgner, P. L., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. *Nature* *337*, 387-388.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nucleic Acids Res.* *15*, 8125-8148.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Mol. Cell. Biol.* *7*, 4125-4129.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* *6*, 742-751.
- Southern, P. J., and Berg, P. (1982). Transformation of Mammalian Cells to Antibiotic Resistance with a Bacterial Gene Under Control of the SV40 Early Region Promoter. *J. Molec. Appl. Gen.* *1*, 327-339.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* *11*, 223-232.

©1998-2003, 2010 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.