

# pVAX1<sup>©</sup>

**Catalog no. V260-20**

**Version C**  
111110  
25-0256



# Table of Contents

General Information .....	v
<b>Methods .....</b>	<b>1</b>
Overview .....	1
Cloning into pVAX1 <sup>®</sup> .....	2
Transient Transfection .....	4
<b>Appendix .....</b>	<b>5</b>
pVAX1 <sup>®</sup> Vector .....	5
pVAX1 <sup>®</sup> //lacZ Map .....	7
Technical Service .....	8
Purchaser Notification .....	10
References .....	11



## General Information

---

**Contents**

20 µg each of pVAX1<sup>®</sup> and pVAX1<sup>®</sup>/*lacZ*, lyophilized in TE, pH 8.0

---

**Shipping/Storage**

Lyophilized plasmids are shipped at room temperature and should be stored at -20°C.

---



# Methods

## Overview

---

### Introduction

pVAX1<sup>®</sup> is a 3.0 kb plasmid vector designed for use in the development of DNA vaccines. The vector was constructed to be consistent with the Food and Drug Administration (FDA) document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Disease Indications”, published December 22, 1996 (see **FDA “Points to Consider”** below). Features of the vector allow high-copy number replication in *E. coli* and high-level transient expression of the protein of interest in most mammalian cells (see page 5-6 for more information). The vector contains the following elements:

- Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells
- Bovine growth hormone (BGH) polyadenylation signal for efficient transcription termination and polyadenylation of mRNA
- Kanamycin resistance gene for selection in *E. coli*

The control plasmid, pVAX1<sup>®</sup>/*lacZ*, is included for use as a positive control for transfection and expression in the cell line of choice.

---

### Experimental Outline

Use the following outline to clone and express your gene of interest in pVAX1<sup>®</sup>.

- Consult the multiple cloning site described on page 3 to design a strategy to clone your gene into pVAX1<sup>®</sup>.
  - Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on LB plates containing 50 µg/ml kanamycin.
  - Analyze your transformants for the presence of insert by restriction digestion.
  - Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in the proper orientation.
  - Transfect your construct into the mammalian cell line of choice and test for transient expression of your protein of interest.
- 

### FDA “Points to Consider”

pVAX1<sup>®</sup> was constructed by modifying the vector, pcDNA<sup>™</sup>3.1, to accommodate the following considerations put forth by the FDA Center for Biologics Evaluation and Research (CBER) in its document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Diseases Indications” (Docket no. 96N-0400).

- Sequences not necessary for replication in *E. coli* or for expression of recombinant protein in mammalian cells were removed to limit DNA sequences with possible homology to the human genome and to minimize the possibility of chromosomal integration.
  - The kanamycin resistance gene was substituted for the ampicillin resistance gene because aminoglycoside antibiotics are less likely to elicit allergic responses in humans.
-

# Cloning into pVAX1<sup>®</sup>

---

## Introduction

A diagram is provided on the next page to help you design a cloning strategy for ligating your gene of interest into pVAX1<sup>®</sup>. General considerations for transformation are listed below.

---

## General Molecular Biology Techniques

For help with DNA ligation, *E. coli* transformation, restriction enzyme analysis, DNA sequencing, and DNA biochemistry, please refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) (See **References**, page 11).

---

## *E. coli* Strain for Transformation

Many *E. coli* strains are suitable for the propagation of this vector, including TOP10, DH5 $\alpha$ , and DH10. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*<sup>1</sup>) and endonuclease A deficient (*endA*).

For your convenience, TOP10 is available as chemically competent or electrocompetent cells from Invitrogen.

Description	Reactions	Efficiency	Catalog no.
One Shot <sup>®</sup> TOP10 Chemically Competent	20 reactions	1 x 10 <sup>9</sup>	C4040-03
One Shot <sup>®</sup> TOP10 Electrocomp <sup>™</sup>	20 reactions	1 x 10 <sup>9</sup>	C4040-50
TOP10 Electrocomp <sup>™</sup>	20 reactions	1 x 10 <sup>9</sup>	C664-55

---

## Transformation Method

You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids.

---

## Maintenance of pVAX1<sup>®</sup>

To propagate and maintain the pVAX1<sup>®</sup> plasmid, resuspend the vector in 20  $\mu$ l sterile water to prepare a 1  $\mu$ g/ $\mu$ l stock solution. Store the stock solution at -20°C.

Use this stock solution to transform a *recA*<sup>1</sup>, *endA* *E. coli* strain like TOP10, DH5 $\alpha$ , or equivalent. Select transformants on LB plates containing 50  $\mu$ g/ml kanamycin. Be sure to prepare a glycerol stock of your plasmid-containing *E. coli* strain for long-term storage (see page 3).

---

## Cloning Considerations

pVAX1<sup>®</sup> is a nonfusion vector. Your insert must contain a Kozak translation initiation sequence and an initiation codon (ATG) for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). An example of a Kozak consensus sequence is provided below. Please note that other sequences are possible (see references above), but the A at position -3 and the G at position +4 are the most critical (shown in bold). The ATG initiation codon is shown underlined.

**ANNATGG**

Your insert must also contain a stop codon for proper termination of your gene. Please note that the *Xba* I site contains an internal stop codon (TCTAGA).

---

*continued on next page*





# Transient Transfection

---

## Introduction

Once you have verified that your gene is cloned in the correct orientation and contains an initiation ATG and a stop codon, then you are ready to transiently transfect your mammalian cell line of choice to check for protein expression. We recommend that you include the positive control vector and a mock transfection (negative control) to evaluate your results.

---

## Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells and salt will interfere with lipids, decreasing transfection efficiency. We recommend isolating plasmid DNA using the S.N.A.P.<sup>™</sup> MiniPrep Kit (10-15 µg DNA, Catalog no. K1900-01), the S.N.A.P.<sup>™</sup> MidiPrep Kit (10-200 µg DNA, Catalog no. K1910-01), or CsCl gradient centrifugation.

---

## Methods of Transfection

For established cell lines (e.g. HeLa), please consult original references or the supplier of your cell line for the optimal method of transfection. It is recommended that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (**Reference** section, page 11).

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1987; Felgner *et al.*, 1989), and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). Invitrogen offers the Calcium Phosphate Transfection Kit for mammalian transfection.

Catalog No.	Description	Quantity
K2780-01	Calcium Phosphate Transfection Kit	75 reactions

---

## Positive Control

pVAX1<sup>®</sup>/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 7). It may be used to optimize transfection conditions for your cell line. The gene encoding β-galactosidase is expressed in mammalian cells under the control of the CMV promoter. A successful transfection will result in β-galactosidase expression that can be easily assayed (see below).

---

## Assay for β-galactosidase Activity

You may assay for β-galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the β-Gal Assay Kit (Catalog no. K1455-01) and the β-Gal Staining Kit (Catalog no. K1465-01) for fast and easy detection of β-galactosidase expression.


---

# Appendix

## pVAX1<sup>®</sup> Vector

### Map of pVAX1<sup>®</sup>

The figure below summarizes the features of the pVAX1<sup>®</sup> vector. The sequence for pVAX1<sup>®</sup> is available for downloading from our World Wide Web site (<http://www.invitrogen.com>) or from Technical Service (see page 8).



**Comments for pVAX1<sup>®</sup>:**  
2999 bp

CMV promoter: bases 137-724  
T7 promoter/priming site: bases 664-683  
Multiple cloning site: bases 696-811  
BGH reverse priming site: bases 823-840  
BGH polyadenylation signal: bases 829-1053  
Kanamycin resistance gene: bases 1226-2020  
pUC origin: bases 2320-2993

*continued on next page*

## pVAX1<sup>®</sup> Vector, continued

### Features of pVAX1<sup>®</sup>

pVAX1<sup>®</sup> (2999 bp) contains the following elements. All features have been functionally tested.

<b>Feature</b>	<b>Benefit</b>
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Permits efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
Multiple cloning site	Allows insertion of your gene and facilitates cloning
BGH reverse priming site	Permits sequencing through the insert
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992)
Kanamycin resistance gene	Selection of vector in <i>E. coli</i> (Davies and Smith, 1978)
pUC origin (pUC-derived)	High-copy number replication and growth in <i>E. coli</i> (Bolivar <i>et al.</i> , 1977; Bolivar <i>et al.</i> , 1977)


# pVAX1<sup>®</sup>//lacZ Map

## Description

pVAX1<sup>®</sup>//lacZ is a 6050 bp control vector containing the gene for β-galactosidase. The vector was constructed by cloning a 3.1 kb *Pst* I-*Xba* I fragment containing the *lacZ* gene into the *Pst* I-*Xba* I site of pVAX1<sup>®</sup>.

## Map of Control Vector

The figure below summarizes the features of the pVAX1<sup>®</sup>//lacZ vector. **The complete nucleotide sequence for pVAX1<sup>®</sup>//lacZ is available for downloading from our World Wide Web site (<http://www.invitrogen.com>) or by contacting Technical Service.** See the next page for more information.



### Comments for pVAX1<sup>®</sup>//lacZ: 6050 bp

CMV promoter: bases 137-724  
T7 promoter/priming site: bases 664-683  
LacZ ORF: bases 773-3829  
BGH reverse priming site: bases 3874-3891  
BGH polyadenylation signal: bases 3880-4104  
Kanamycin resistance gene: bases 4277-5071  
pUC origin: bases 5371-6044

# Technical Service

---

## World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

## Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page ([www.invitrogen.com](http://www.invitrogen.com)).

### Corporate Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
E-mail:  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

### Japanese Headquarters:

Invitrogen Japan K.K.  
Nihonbashi Hama-Cho Park Bldg. 4F  
2-35-4, Hama-Cho, Nihonbashi  
Tel: 81 3 3663 7972  
Fax: 81 3 3663 8242  
E-mail: [jpinfo@invitrogen.com](mailto:jpinfo@invitrogen.com)

### European Headquarters:

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK  
Tel: +44 (0) 141 814 6100  
Tel (Toll Free): 0800 5345 5345  
Fax: +44 (0) 141 814 6117  
E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

## MSDS Requests

To request an MSDS, visit our Web site ([www.invitrogen.com](http://www.invitrogen.com)) and follow the instructions below.

1. On the home page, go to the left-hand column under 'Technical Resources' and select 'MSDS Requests'.
  2. Follow instructions on the page and fill out all the required fields.
  3. To request additional MSDSs, click the 'Add Another' button.
  4. All requests will be faxed unless another method is selected.
  5. When you are finished entering information, click the 'Submit' button. Your MSDS will be sent within 24 hours.
- 

*Continued on next page*

## Technical Service, continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

## Purchaser Notification

---

**Limited Use Label  
License No. 358:  
Research Use  
Only**

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

---



## References

---

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Bolivar, F., Rodriguez, R. L., Betlach, M. C., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. I. Ampicillin-Resistant Derivatives of the Plasmid pMB9. *Gene* *2*, 75-93.
- Bolivar, F., Rodriguez, R. L., Greene, P. J., Betlach, M. C., Heyneker, H. L., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. II. A Multipurpose Cloning System. *Gene* *2*, 95-113.
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nuc. Acids Res.* *15*, 1311-1326.
- Davies, J., and Smith, D. I. (1978). Plasmid-Determined Resistance to Antimicrobial Agents. *Ann. Rev. Microbiol.* *32*, 469-518.
- Felgner, P. L., Gadek, T. R., Holm, M., Roman, R., Chan, H. W., Wenz, M., Northrop, J. P., Ringold, G. M., and Danielsen, M. (1987). Lipofectin: A Highly Efficient, Lipid-mediated DNA-transfection Procedure. *Proc. Natl. Acad. Sci. USA* *84*, 7413-7417.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nuc. Acids Res.* *15*, 8125-8148.
- Kozak, M. (1991). An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biol.* *115*, 887-903.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Mol. Cell. Biol.* *7*, 4125-4129.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* *6*, 742-751.