

pYD1 Yeast Display Vector Kit

Catalog no. V835-01

Version D

December 10, 2002

25-0259



www.invitrogen.com
tech_service@invitrogen.com

Table of Contents

Table of Contents.....	iii
Important Information	v
Purchaser Notification	vi
Introduction	1
Overview	1
Methods	3
Cloning into pYD1	3
Yeast Transformation	7
Expression and Display of Protein Fusion.....	9
Appendix.....	13
Recipes.....	13
pYD1 Vector	15
Small-Scale Yeast Transformation	17
Technical Service	18
References	20

Important Information

Contents

20 µg of pYD1, lyophilized in TE, pH 8
EBY100 glycerol stock, 0.5 ml
2 µg pYD1 Forward sequencing primer, lyophilized
2 µg pYD1 Reverse sequencing primer, lyophilized

Shipping/Storage

The pYD1 Vector Kit is shipped on dry ice. Upon receipt:

- Store the glycerol stock at -80°C
 - Store the lyophilized plasmid and primers at -20°C
-

Information for European Customers

The EGY100 yeast strain is genetically modified and contains the plasmid, pIU211 stably integrated into the genome (Boder and Wittrup, 1997). As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

Purchaser Notification

**Limited Use Label
License No. 88:
Yeast Display
Expression Kit**

The Yeast Display Expression Kit is covered by filed US Patent Applications and is sold under license from The Board of Trustees of the University of Illinois for life science research only; not for any diagnostic, therapeutic or commercial purposes (other than (d) and (e) below) for evaluation for one (1) year following the date of purchase of the Expression Kit. Neither the Expression Kit, nor the materials contained therein, may be used for any commercial purpose without a license for such purposes from The Board of Trustees of the University of Illinois, Urbana, Illinois. Commercial purposes include: (a) any use of Expression Products in a Commercial Product; (b) any use of Expression Products in the manufacture of a Commercial Product; (c) any sale of Expression Products; (d) any use of Expression Products or the Expression Kit to facilitate or advance any research or development of a Commercial Product; and (e) any use of Expression Products or the Expression Kit to facilitate or advance any research or development program the results of which will be applied to the development of Commercial Products. "Expression Product(s)" means products expressed with or through use of the Expression Kit, or with or through the use of any vectors or host strains in the Expression Kit. "Commercial Product(s)" means any product intended for sale or commercial use. Commercial entities may use this Expression Kit for the commercial purposes specified in (d) and (e) above solely for the evaluation of possible commercial applications and only for a period of one (1) year from the date of purchase of the Expression Kit, during which time they must contact the University of Illinois regarding a commercial license for use of the Expression Kit for commercial purposes, or thereafter immediately terminate their use of, and destroy any remaining components of, the Expression Kit, and so notify the University of Illinois in writing. Access to the Expression Kit must be limited solely to those officers, employees and students of your institution who need access to perform the above-described research or evaluation. You may not distribute the Expression Kit or the vectors or host strains contained in it, nor modifications or alterations thereof, to others, even those within your own institution, without prior written approval from the University of Illinois. Failure to comply with these terms and conditions of use shall result in immediate termination of your limited right to use the Expression Kit at which point you must destroy all Expression Kits in your control, and so notify Invitrogen Corporation in writing. For information regarding commercial licenses, contact: University of Illinois, Research and Technology Management Office, Fourth Floor, Swanlund Building, 601 East John Street, Champaign, IL 61820. Tel: (217) 333-7862, Fax: (217) 244-3716.

**Limited Use Label
License No. 22:
Vectors and
Clones Encoding
Histidine Hexamer**

This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

**Limited Use Label
License No. 141:
Expression of
Polypeptides in
Yeast**

This product is the subject of one or more of U.S. Patent Nos. 5,618,676, 5,854,018, 5,856,013, 5,919,651 and foreign equivalents. Rights to use this product are limited to academic research use only. For-profit entities are required to obtain a separate license from Washington Research Foundation to utilize this product for any use. Washington Research Foundation, 2815 Eastlake Avenue East, Suite 300, Seattle, Washington 98102. Tel: 206-336-5600. Fax: 206-336-5615.

Introduction

Overview


Introduction

pYD1 is a 5.0 kb expression vector designed for expression, secretion, and display of proteins on the extracellular surface of *Saccharomyces cerevisiae* cells. Features of this vector allow regulated expression, secretion, and detection of expressed proteins (see pages 15-16 for more information). The vector contains the following elements:

- *AGA2* gene from *Saccharomyces cerevisiae*. This gene encodes one of the subunits of the **a**-agglutinin receptor. Fusion of the gene of interest to *AGA2* allows secretion and display of the protein of interest.
- *GAL1* promoter for regulated expression of the *AGA2* gene fusion.
- Xpress™ epitope and V5 epitope for detection of the displayed protein.
- Polyhistidine (6xHis) tag for detection and possible purification on metal chelating resin.
- *TRP1* gene for selection in *Saccharomyces cerevisiae*.
- *CEN6/ARS4* for stable, episomal replication in yeast.
- Ampicillin resistance gene and the pUC origin for selection and replication in *E. coli*.

Agglutinin Receptor

The native **a** and **α** agglutinin receptors are believed to act as adhesion molecules to stabilize cell-cell interactions during mating and to facilitate fusion between the **a** and **α** haploid yeast cells. The Yeast Display System uses the **a**-agglutinin receptor of *S. cerevisiae* to display foreign proteins on the cell surface (Boder and Wittrup, 1997). The **a**-agglutinin receptor consists of two subunits encoded by the *AGA1* and *AGA2* genes (see diagram below). The Aga1 protein (Aga1p, 725 amino acids) is secreted from the cell and becomes covalently attached to β-glucan in the extracellular matrix of the yeast cell wall (Lu *et al.*, 1995). The Aga2 protein (Aga2p, 69 amino acids) binds to Aga1p through two disulfide bonds and after secretion remains attached to the cell through its contact with Aga1p (Cappellaro *et al.*, 1994; Cappellaro *et al.*, 1991; Roy *et al.*, 1991). The N-terminal portion of Aga2p is required for attachment to Aga1p, while proteins and peptides can be fused to the C-terminus for presentation on the yeast cell surface (Boder and Wittrup, 1997).



continued on next page

Overview, continued

Experimental Outline

Use the following outline to clone and express your gene of interest in pYD1

- Consult the multiple cloning site on page 5 to clone your gene in frame with the N-terminal leader peptide encoding Aga2p, a Gly-Ser linker, and the Xpress™ epitope and the C-terminal peptide encoding the V5 epitope and the 6xHis tag (optional).
 - Ligate your insert and transform into *E. coli*. Select transformants on LB containing 50-100 µg/ml ampicillin.
 - Analyze your transformants for the presence of insert by restriction digestion.
 - Select a transformant with the correct restriction pattern and sequence to confirm that your gene is cloned in frame with the N-terminal and C-terminal peptides.
 - Transform your construct into the yeast strain EBY100 and select transformants on minimal medium without tryptophan.
 - Select a transformant and grow overnight in glucose minimal medium.
 - Induce expression with galactose medium and test for display of your protein.
-

Applications

Yeast display vectors offer an alternative to other display technologies such as phage display, bacterial surface display, and the yeast two-hybrid method. Recently, a yeast display vector has been used to apply library methods for *in vitro* antibody affinity maturation (Boder and Wittrup, 1998; Boder and Wittrup, 1997; Kieke *et al.*, 1997). A yeast display system is particularly useful for mammalian cell surface and secreted proteins (e.g., receptors, cytokines) that require endoplasmic reticulum-specific posttranslational processing for efficient folding and activity. In addition to antibody affinity maturation, other potential applications include:

- Isolation of binding domains from cDNA expression libraries
 - Isolation of mutant receptors or ligands
-

Methods

Cloning into pYD1

Introduction

A diagram is provided on page 5 to help you ligate your gene of interest in frame with the N-terminal and the C-terminal peptide (optional). General considerations for cloning and transformation are listed below.

General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, see *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

E. coli Strain

Many *E. coli* strains are suitable for the growth of this vector including TOP10 and DH5 α . We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*) and endonuclease A deficient (*endA*).

For your convenience, TOP10 is available as chemically competent or electrocompetent cells from Invitrogen.

Item	Quantity	Catalog no.
One Shot™ TOP10 Chemically Competent <i>E. coli</i>	21 x 50 μ l	C4040-03
Electrocomp™ TOP10	5 x 80 μ l	C664-55

Transformation Method

Use your method of choice to transform pYD1 containing your insert into *E. coli*. Chemical transformation is the most convenient method for most researchers. Electroporation is the most efficient and the method of choice for large plasmids (>10 kb).

Maintenance of pYD1

In order to propagate and maintain pYD1, we recommend that you transform the plasmids into *E. coli* and prepare glycerol stocks for long-term storage. To transform plasmids into *E. coli*:

- Resuspend the lyophilized vector in 20 μ l sterile water to make a 1 μ g/ μ l stock solution. Store the stock solution at -20°C.
- Use this stock solution (10-100 ng) to transform a *recA*, *endA* *E. coli* strain like TOP10 or DH5 α .
- Select transformants on LB plates containing 50-100 μ g/ml ampicillin.

To prepare glycerol stocks, see page 6.

continued on next page

Cloning into pYD1, continued

Detection of Recombinant Proteins

Before you clone your gene of interest into pYD1, it is best to consider how you will detect your recombinant fusion on the surface of yeast cells. You may detect your protein by activity assay or by fluorescent staining. Usually fluorescent staining with the appropriate antibody is used. If you do not have an antibody to your protein, there are several antibodies from Invitrogen which you can utilize. Ordering information is provided below.

Antibody	Epitope	Catalog no.
Anti-Xpress™ Antibody	Detects 8 amino acid epitope	R910-25
Anti-Xpress™-HRP Antibody	DLYDDDDK	R911-25
Anti-Xpress™-FITC Antibody		R913-25
Anti-V5 Antibody	Detects 14 amino acid epitope derived from SV5 (Southern <i>et al.</i> , 1991)	R960-25
Anti-V5-HRP Antibody		R961-25
Anti-V5-AP Antibody	GKPIPPLLGLDST	R962-25
Anti-V5-FITC Antibody		R963-25
Anti-His(C-term) Antibody	Detects the C-terminal polyhistidine tag (requires free carboxyl group for detection) (Lindner <i>et al.</i> , 1997)	R930-25
Anti-His(C-term)-HRP Antibody		R931-25
Anti-His(C-term)-AP Antibody	HHHHHH-COOH	R932-25
Anti-His(C-term)-FITC Antibody		R933-25

Cloning into pYD1

pYD1 is a fusion vector requiring that you clone your gene of interest in frame with the *AGA2* gene and the C-terminal V5 epitope/polyhistidine tag (optional). For proper expression, first determine if any restriction sites are appropriate to preserve the reading frame at **BOTH** the 5' and the 3' ends. It may be necessary to PCR your gene product to create a fragment with the appropriate restriction sites to clone in frame at both ends. Carefully inspect your gene and the multiple cloning site before cloning your gene of interest. See the next page for details of the multiple cloning site.



Note

If you wish to express and display your protein without the C-terminal tag, include a termination codon in your gene of interest. Note that if you include a termination codon, you will be unable to detect the fusion protein with the Anti-V5 Antibody or Anti-His (C-term) Antibody.

Other Cloning Strategies

If you wish to clone your gene closer to the *AGA2* gene, the *Hind* III or the *Pst* I sites between the *AGA2* gene and the DNA encoding the Gly-Ser linker may be used. There is also an *Nhe* I site between the Gly-Ser linker and the Xpress™ epitope which can be used. Several different cloning strategies may be tested to ensure that your protein is displayed properly for your particular application.

continued on next page

Cloning into pYD1, continued

Multiple Cloning Site of pYD1

Below is the multiple cloning site for pYD1. Restriction sites are labeled to indicate the cleavage site. The multiple cloning site has been confirmed by sequencing and functional testing. **The complete sequence is available by downloading it from our Web site (www.invitrogen.com) or by calling Technical Service (page 18).**

```

transcriptional start ----->
320 AAACATGCAT AACCACTTTA ACTAATACTT TCAACATTTT CGGTTTGTAT TACTTCTTAT TCAAATGTAA TAAAAGTATC

400 AACAAAAAAT TGTTAATATA CCTCTATACT TTAACGTCAA GGAGAAAAAA CCCCGGATCG GACTACTAGC AGCTGTAATA

T7 promoter/priming site -----
480 CGACTCACTA TAGGGAATAT TAAGCTAATT CTACTTCATA CATTTTCAAT TAAG ATG CAG TTA CTT CGC TGT TTT
Met Gln Leu Leu Arg Cys Phe

Start of Aga2 mature peptide -----
565 TCA ATA TTT TCT GTT ATT GCT TCA GTT TTA GCA CAG GAA CTG ACA ACT ATA TGC GAG CAA ATC CCC
Ser Ile Phe Ser Val Ile Ala Ser Val Leu Ala Gln Glu Leu Thr Thr Ile Cys Glu Gln Ile Pro
Cleavage site ^

621 TCA CCA ACT TTA GAA TCG ACG CCG TAC TCT TTG TCA ACG ACT ACT ATT TTG GCC AAC GGG AAG GCA
Ser Pro Thr Leu Glu Ser Thr Pro Tyr Ser Leu Ser Thr Thr Thr Ile Leu Ala Asn Gly Lys Ala

pYD1 Forward priming site -----
687 ATG CAA GGA GTT TTT GAA TAT TAC AAA TCA GTA ACG TTT GTC AGT AAT TGC GGT TCT CAC CCC TCA
Met Gln Gly Val Phe Glu Tyr Tyr Lys Ser Val Thr Phe Val Ser Asn Cys Gly Ser His Pro Ser

End of Aga2 mature peptide ----- Hind III Pst I*
753 ACA ACT AGC AAA GGC AGC CCC ATA AAC ACA CAG TAT GTT TTT AAG CTT CTG CAG GCT AGT GGT GGT
Thr Thr Ser Lys Gly Ser Pro Ile Asn Thr Gln Tyr Val Phe Lys Leu Leu Gln Ala Ser Gly Gly

Gly-Ser Linker ----- Nhe I
819 GGT GGT TCT GGT GGT GGT GGT TCT GGT GGT GGT GGT TCT GCT AGC ATG ACT GGT GGA CAG CAA ATG
Gly Gly Ser Gly Gly Gly Gly Ser Gly Gly Gly Ser Ala Ser Met Thr Gly Gly Gln Gln Met

Xpress™ epitope ----- Acc65 I Kpn I BamH I BstX I EcoR I Pst I*
885 GGT CGG GAT CTG TAC GAC GAT GAC GAT AAG GTA CCA GGA TCC AGT GTG GTG GAA TTC TGC AGA TAT
Gly Arg Asp Leu Tyr Asp Asp Asp Asp Lys Val Pro Gly Ser Ser Val Val Glu Phe Cys Arg Tyr
Enterokinase recognition site ^ Enterokinase cleavage site

BstX I Not I Xho I Xba I PspOM I Apa I BstB I ----- V5 epitope
951 CCA GCA CAG TGG CGG CCG CTC GAG TCT AGA GGG CCC TTC GAA GGT AAG CCT ATC CCT AAC CCT CTC
Pro Ala Gln Trp Arg Pro Leu Glu Ser Arg Gly Pro Phe Glu Gly Lys Pro Ile Pro Asn Pro Leu

Polyhistidine region ----- Pme I
1017 CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC CAT CAC CAT TGA GTTTA AACCCGCTGA
Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His His His His ***

pYD1 Reverse priming site -----
1080 TCTGATAACA ACAGTGTA GAAGTAAATA TCGACTTTGT TCCCACTGTA CTTTTAGCTC GTACAAAATA CAATATACTT

1160 TTCATTCTC CGTAAACAAC ATGTTTCC C ATGTAATATC CTTTTCTATT TTTCGTTCCG TTACCAACTT TACACATACT

```

*The two Pst I sites are unique to the multiple cloning site.

continued on next page

Cloning into pYD1, continued

E. coli Transformation

Transform your ligation mixtures into a competent *recA*, *endA* *E. coli* strain and select on LB plates containing 50-100 µg/ml ampicillin. Select 10 colonies and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the pYD1 Forward and Reverse primers included in the kit to confirm that your gene is in frame with the N-terminal peptide and the C-terminal peptide.

Sequence of Primers

The table below provides the sequences and number of pmoles supplied for the pYD1 Forward and pYD1 Reverse sequencing primers. If you need to purchase more primer, ordering information is also provided.

Primer	Sequence	Pmoles	Catalog no.
pYD1 Forward	5'-AGTAACGTTTGTTCAGTAATTGC-3'	297	N842-02
pYD1 Reverse	5'-GTCGATTTTGTTCATCTACAC-3'	301	N843-02

Preparing a Glycerol Stock

Once you have identified the correct clone, be sure to purify the colony and make a glycerol stock for long-term storage. It is also a good idea to keep a DNA stock of your plasmid at -20°C as well.

- Streak the original colony out on an LB plate containing 50-100 µg/ml ampicillin. Incubate the plate at 37°C overnight.
- Isolate a single colony and inoculate into 1-2 ml of LB containing 50-100 µg/ml ampicillin.
- Grow the culture to mid-log phase ($OD_{600} = 0.5-0.7$).
- Mix 0.85 ml of culture with 0.15 ml of sterile glycerol and transfer to a cryovial.
- Store at -80°C.

Yeast Transformation

Introduction

Once you have confirmed that your construct is in the correct orientation and fused in frame with the N-terminal peptide and the C-terminal peptide (if desired), you are ready to transform the yeast strain EBY100.

Basic Yeast Molecular Biology

The user should be familiar with basic yeast molecular biology and microbiological techniques. Refer to *Current Protocols in Molecular Biology*, Unit 13 (Ausubel *et al.*, 1994) and the *Guide to Yeast Genetics and Molecular Biology* (Guthrie and Fink, 1991) for information on preparing yeast media and handling yeast.

Positive Control

We recommend that you transform the empty vector (pYD1) into EBY100 for use as a positive control. pYD1 expresses Aga2p fused to the Xpress™, V5, and polyhistidine epitopes. Expression and detection of this fusion protein by fluorescent staining will demonstrate that the system is working.

EBY100

EBY100 expresses the *AGAI* gene under control of the *GALI* promoter. It was created by integrating the vector pIU211 into the *AGAI* locus of the *Saccharomyces cerevisiae* strain BJ5465 (*MATa ura 3-52 trp 1 leu2Δ1 his3Δ200 pep4:HIS3 prb1Δ1.6R can1 GAL*) (Boder and Wittrop, 1997). pIU211 contains the *AGAI* gene regulated by the *GAL* promoter and a *URA3* selectable marker. The vector was linearized for integration at *AGAI* using *Bst*W I. *Bst*W I cuts in the *AGAI* gene, ensuring that the construct integrates at the *AGAI* locus.

EBY100 (Catalog no. C839-00) is available separately.

Transformation Method

You may use your own method of choice to transform EBY100.

If you do not have a method for transformation, we recommend the *S.c.* EasyComp™ Kit (Catalog no. K5050-01). This kit rapidly prepares competent *Saccharomyces cerevisiae* in less than 30 minutes. In addition, the cells can be frozen for future use.

Alternatively, we have provided a small-scale transformation protocol utilizing lithium acetate in the Appendix (page 17).

Before Starting

To transform EBY100 you will need the following reagents. For media recipes, refer to page 13.

- Minimal dextrose plates containing leucine and tryptophan
 - YPD medium
 - Minimal dextrose plates containing leucine
 - 20% Glucose
 - 0.5-5 µg plasmid DNA (DNA may be isolated from *E. coli* using your method of choice)
-

continued on next page

Yeast Transformation, continued

Transformation of EBY100

Follow the general steps below to transform EBY100.

1. Take the glycerol stock of EBY100 and streak out on a minimal dextrose plate containing leucine and tryptophan. Incubate at 30°C until colonies appear (1-2 days).
2. Use your method of choice to prepare competent cells and transform with pYD1 or pYD1 containing your gene of interest.
3. Plate the transformation reaction (50-150 μ l) on minimal dextrose plates containing leucine. Incubate the plates at 30°C for 2 to 4 days until single colonies appear.

Once you have transformants, you are ready to test for display of your fusion protein. Proceed to the next page to induce expression of your fusion protein.

Expression and Display of Protein Fusion

Introduction

To express your fusion protein in EBY100, transformants are grown in glucose-containing medium overnight, then switched to medium containing galactose for induction of expression. Cells are grown in galactose until display of the desired protein fusion is detected by antibody staining. This is usually a function of the protein fusion and is determined empirically. Follow the instructions below to express and detect your protein fusion.

Controls

We recommend including the following controls to evaluate your experiment.

Control	Type of Control	Purpose
EBY100/pYD1	Positive	Expression and detection of the Aga2 fusion protein demonstrates that the system is working properly.
EBY100	Negative	Controls for background resulting from the cells
EBY100/pYD1 + gene of interest (uninduced)	Negative	Controls for leaky expression. (This will be the zero time point for the induction time course on page 10)

Before Starting

Be sure to have following reagents and equipment on hand. Refer to page 13 for media recipes.

- Shaking incubator (adjustable from 20°C to 30°C) **Note:** A shaking incubator with a cooling apparatus may be required to maintain 20°C consistently
 - 50 ml conical centrifuge tube or 25-50 ml culture flask
 - Spectrophotometer
 - YNB-CAA medium
 - 20% glucose
 - 20% galactose
 - 1X Phosphate Buffered Saline (PBS)
 - 1X PBS, 1 mg/ml BSA
 - Anti-mouse IgG conjugated with FITC (if you are using a polyclonal antibody, use anti-rabbit IgG conjugated with FITC) **Note:** Other fluorophores may be suitable.
-

continued on next page

Expression and Display of Protein Fusion, continued

Induction

Grow untransformed EBY100, EBY100/pYD1, and EBY100/pYD1 containing your gene of interest as follows.

1. Inoculate a single yeast colony into 10 ml YNB-CAA containing 2% glucose and grow overnight at 30°C with shaking.
2. Read the absorbance of the cell culture at 600 nm. The OD₆₀₀ should be between 2 and 5. If the OD₆₀₀ is below 2 or over 5, refer to the table below.

If the OD ₆₀₀ is....	Then you may...
below 2,	either continue growing the cells until the OD ₆₀₀ reaches 2 or proceed to Step 3.
over 5,	proceed to Step 3.

Note: In general, OD₆₀₀ readings less than 2 will decrease the number of cells displaying fusion proteins. OD₆₀₀ readings greater than 5 will delay induction of expression of the displayed protein.

3. Centrifuge the cell culture at 3000-5000 x g for 5-10 minutes at room temperature.
4. Resuspend the cell pellet in YNB-CAA medium containing 2% galactose to an OD₆₀₀ of 0.5 to 1. This is to ensure that the cells continue to grow in log-phase. For example, if the OD₆₀₀ is 2 from Step 2, resuspend the cells in 20 to 40 ml of medium.
5. Immediately remove a volume of cells equivalent to 2 OD₆₀₀ units. For an OD₆₀₀ of 0.5, remove 4 ml and place on ice. This is your zero time point.
6. Incubate the cell culture at 20-25°C with shaking. **Note:** In general, more cells will display the protein fusion at 20°C. However, if you do not have an incubator that maintains 20°C, try 25°C.
7. Assay the cell culture over a 48-hour time period (i.e. 0, 12, 24, 36, 48 hours) to determine the optimal induction time for maximum display. For each time point, read the OD₆₀₀ and remove a volume of cells that is equivalent to 2 OD₆₀₀ units (see Step 5, above). Proceed to **Staining of Displayed Proteins**, below.

Staining of Displayed Proteins

For each time point, assay untransformed EBY100, EBY100/pYD1, and EBY100/pYD1 containing your gene of interest. Time points may be processed as they are collected or placed on ice and stored at +4°C until all time points are collected. **Do not freeze cells.**

1. Take your time points from Steps 4 and 6, above, and centrifuge at 3000-5000 x g for 5-10 minutes at +4°C.
2. Resuspend the cells in 1X PBS and centrifuge as in Step 1.
3. Remove the PBS and resuspend the cell pellet in 250 µl of 1X PBS, 1 mg/ml BSA, and 1 µg antibody.
4. Incubate on ice for 30 minutes with occasional mixing.
5. Centrifuge the cells at 3000-5000 x g for 5-10 minutes at +4°C.
6. Wash the cells once with 1X PBS.
7. Resuspend the cells in 1X PBS, 1 mg/ml BSA, and 1 µg anti-mouse IgG conjugated with fluorescein isothiocyanate (FITC) or any other suitable fluorophore.

continued on next page

Expression and Display of Protein Fusion, continued

Staining of Displayed Proteins, continued

8. Incubate on ice for 30 minutes in the dark and invert the tube occasionally.
 9. Wash the cells two times with 1 ml of 1X PBS (no BSA).
 10. Resuspend the cells in 40 μ l 1X PBS and spot 5-10 μ l on a slide and add a coverslip for microscopic observation under UV light. Record your observations or take a picture.
-

Evaluating Your Results

The Aga2p fusion can be detected as early as 4 hours after the switch to galactose and is still detectable after 72 hours in galactose medium. Optimal detection occurs between 12 and 48 hours.

You should see little or no fluorescence with untransformed EBY100, uninduced EBY100/pYD1 (zero time point), or uninduced EBY100/pYD1 containing your gene of interest (zero time point).

If you have trouble expressing and displaying your protein, see the next page.

Examples of Displayed Proteins

Growth conditions (i.e. density of cells, temperature) and the type of protein are the major variables affecting expression and display. Displayed proteins can be detected as early as 4 hours and as late as 72 hours after the shift from glucose to galactose. Detection of an Aga2-GFP fusion using the V5 antibody is summarized below.

Time (hours)	Fluorescence
2	Barely detectable above background
4	Detectable, but faint
8	Much better
20-30	Optimal fluorescence

Other proteins, such as an Aga2-T-cell receptor fusion, required 36-48 hours of galactose induction for optimal display (K. D. Wittrup, personal communication).

Percentage of Cells Displaying Protein

Under the conditions described above, 10% to 50% of the cells are easily detectable as displaying the fusion protein. If Fluorescence Activated Cell Sorting (FACS) is used to detect fluorescence, up to 70% of the total cells can be detected as displaying the fusion protein in question (Boder and Wittrup, 1997). Note that this is because of the increased sensitivity of FACS, and the fact that induction of expression was done at 20°C rather than 25°C.

Protein Purification

At this time, we have not tried to purify displayed proteins utilizing the polyhistidine tag and metal-chelating resin. Displayed proteins tend to be glycosylated which may or may not interfere with purification. It is possible to release the Aga2p protein fusion from the cell wall by treatment with 100 μ M dithiothreitol (DTT) in buffer for two hours. DTT reduces the disulfide bond holding Aga2p to Aga1p. It is possible that this could be exploited to purify the protein of interest.

continued on next page

Expression and Display of Protein Fusion, continued

Troubleshooting

Use the table below to troubleshoot your experiment.

Problem	Reason	Solution
No display detected with construct or positive control	Induction not performed properly	Review the protocol on page 10
	Cells not growing	Test for growth on minimal medium without uracil or tryptophan. EBY100 transformed with pYD1 should be able to grow without uracil or tryptophan. Review media recipes and re-make if necessary.
No display detected with construct, but display is detected with the positive control	Gene is not fused correctly to the <i>AGA2</i> gene	Sequence construct to confirm that gene is cloned in frame with <i>AGA2</i> . Re-make construct if necessary.
Construct expresses and displays protein fusion, but percentage of cells is low	Cells not grown to a high enough density before induction	Grow culture on glucose to a density > 2 before the shift to galactose
	Induction temperature is too high	Try inducing your construct at 20°C. Incubation at higher temperatures can lead to poor surface expression. This may be because of inefficient protein folding and retention in the endoplasmic reticulum (ER)
Construct expresses and displays protein fusion, but it takes over 48 hours	Cell density was greater than 5 prior to the shift to galactose	Do not grow cells as long on glucose

Appendix

Recipes

LB (Luria-Bertani) Medium and Plates

Composition:

1.0% Tryptone
0.5% Yeast Extract
1.0% NaCl
pH 7.0

1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 ml deionized water.
2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
3. Autoclave on liquid cycle for 20 minutes. Allow solution to cool to 55°C and add antibiotic if needed.
4. Store at room temperature or at +4°C.

LB agar plates

1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.
 2. Autoclave on liquid cycle for 20 minutes at 15 psi.
 3. After autoclaving, cool to ~55°C, add antibiotic (50-100 µg/ml of ampicillin), and pour into 10 cm plates.
 4. Let harden, then invert and store at +4°C, in the dark. Plates are stable for about 2 weeks with antibiotic.
-

YPD

Yeast Extract Peptone Dextrose Medium (1 liter)

1% yeast extract
2% peptone
2% dextrose (D-glucose)

1. Dissolve the following in 900 ml of water:
 - 10 g yeast extract
 - 20 g of peptone
2. Optional: Add 20 g agar, if making plates.
3. Autoclave for 20 minutes on liquid cycle and cool to ~55-60°C.
4. Add 100 ml of 20% dextrose (filter-sterilized).

Store medium at room temperature. The shelf life is approximately one to two months.

continued on next page

Recipes, continued

Minimal Dextrose Plates (with Amino Acids)

0.67% YNB (with ammonium sulfate, **without** amino acids)
2% glucose
0.01% leucine
0.01% tryptophan (optional)
1.5% agar

Note: Plates without tryptophan are required for selection of pYD1 transformants.

1. For 1 liter, dissolve 6.7 g YNB and 15 g agar in 900 ml deionized water.
 2. Autoclave on liquid cycle for 20 minutes.
 3. After autoclaving, allow the medium to cool to ~55-60°C and add 100 ml of 20% glucose (filter-sterilized).
 4. Add 10 ml of 10 mg/ml leucine (filter-sterilized) and/or 10 ml of 10 mg/ml tryptophan (filter-sterilized). Alternatively, plates can be made ahead of time and amino acids spread onto the plate prior to use. Spread 100 µl of a 10 mg/ml amino acid solution onto the plate and let the solution soak into the plate before plating cells.
 5. Store at +4°C. Plates are stable for 1-2 months.
-

YNB-CAA Medium

0.67% YNB (with ammonium sulfate, **without** amino acids)
0.5% Casamino acids
2% glucose or galactose

1. For 1 liter, dissolve 6.7 g YNB and 5 g Casamino acids in 900 ml deionized water.
2. Autoclave on liquid cycle for 20 minutes. **Note:** Casamino acids contain trace amounts of tryptophan that is destroyed upon autoclaving. Lack of tryptophan helps maintain pYD1 during induction.
3. After autoclaving, cool to 55-60°C and add 100 ml of 20% glucose or galactose.
4. Store medium at room temperature or +4°C. Medium is stable for about 1 month.

Alternatively, you may use prepared drop-out medium (i.e. Bio101 HSM-trp-ura) instead of Casamino acids. Omitting uracil helps to maintain the *AGAI* cassette integrated in EBY100.

1X PBS

137 mM NaCl
2.7 mM KCl
10 mM Na₂HPO₄
1.8 mM KH₂PO₄
pH 7.4

1. Dissolve:
 - 8 g NaCl
 - 0.2 g KCl
 - 1.44 g Na₂HPO₄
 - 0.24 g KH₂PO₄in 800 ml deionized water.
 2. Adjust pH to 7.4 with concentrated HCl.
 3. Bring the volume to 1 liter. You may autoclave the solution to increase shelf life.
-

pYD1 Vector

Features of pYD1 pYD1 (5009 bp) contains the following elements. All features have been functionally tested.


Feature	Benefit
<i>GAL1</i> promoter	Permits regulated expression of your recombinant protein (Giniger <i>et al.</i> , 1985; West <i>et al.</i> , 1984)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation
<i>AGA2</i> gene	Allows secretion of the fusion protein
pYD1 Forward priming site	Sequencing of insert
Gly-Ser Linker	Permits spacing of protein of interest from Aga2p
Xpress™ epitope (Asp-Leu-Tyr-Asp-Asp-Asp-Lys)	Allows detection of the fusion protein with the Anti-Xpress™ Antibodies
Multiple cloning site	Allows insertion of your gene of interest
V5 epitope	Allows detection of the fusion protein with the Anti-V5 Antibodies
Polyhistidine epitope (6xHis tag)	Permits detection of pYD1 fusion proteins with the Anti-His(C-term) Antibodies (Lindner <i>et al.</i> , 1997) In addition, it may allow purification of pYD1 fusion proteins by chromatography on metal-chelating resins such as ProBond™
pYD1 Reverse priming site	Permits sequencing through the insert
<i>MATα</i> transcription termination region	Efficient transcription termination and polyadenylation of mRNA
<i>TRP1</i> gene	Allows selection of yeast transformants
<i>CEN6/ARS4</i>	Allows stable, episomal replication and partitioning into daughter cells in yeast
Ampicillin resistance gene (β-lactamase)	Selection in <i>E. coli</i>
pUC origin	High-copy number replication and growth in <i>E. coli</i>

continued on next page

pYD1 Vector, continued

Map of pYD1

The figure below shows the location of features in pYD1. The complete nucleotide sequence for pYD1 is available for downloading from our World Wide Web site (<http://www.invitrogen.com>) or from Technical Service (page 18).



Comments for pYD1 5009 nucleotides

GAL1 promoter: bases 1-451

T7 promoter priming site: bases 475-494

AGA2 ORF: bases 534-794

pYD1 Forward priming site: bases 751-770

Gly/Ser Linker: bases 807-857

Xpress™ epitope: bases 891-914

Multiple cloning site: bases 914-992

V5 epitope: bases 993-1034

Polyhistidine region: bases 1044-1061

pYD1 Reverse priming site: bases 1114-1093 (C)

MAT α transcription termination region: bases 1385-1098 (C)

TRP1 ORF: bases 2107-1433 (C)

CEN6/ARS4: bases 2797-2282 (C)

Ampicillin resistance gene (ORF): bases 2931-3791

pUC origin: bases 3936-4609

Small-Scale Yeast Transformation

Introduction

A small-scale yeast transformation protocol is provided below.

Materials Needed

Be sure to have the following reagents on hand before starting.

- YPD
 - 1X TE
 - 1X LiAc/0.5X TE
 - Denatured sheared salmon sperm DNA
 - Plasmid DNA to be transformed
 - 1X LiAc/40% PEG-3350/1X TE
 - DMSO
 - Selective plates
-

Protocol

1. Inoculate 10 ml of YPD with a yeast colony and shake overnight at 30°C.
 2. Determine the OD₆₀₀ of your overnight culture. Dilute culture to an OD₆₀₀ of 0.4 in 50 ml of YPD and grow an additional 2-4 hours.
 3. Pellet the cells at 2500 rpm and resuspend the pellet in 40 ml 1X TE.
 4. Pellet the cells at 2500 rpm and resuspend pellet in 2 ml of 1X LiAc/0.5X TE.
 5. Incubate the cells at room temperature for 10 minutes.
 6. For each transformation, mix together 1 µg plasmid DNA and 100 µg denatured sheared salmon sperm DNA with 100 µl of the yeast suspension from Step 5.
 7. Add 700 µl of 1X LiAc/40% PEG-3350/1X TE and mix well.
 8. Incubate solution at 30°C for 30 minutes.
 9. Add 88 µl DMSO, mix well, and heat shock at 42°C for 7 minutes.
 10. Centrifuge in a microcentrifuge for 10 seconds and remove supernatant.
 11. Resuspend the cell pellet in 1 ml 1X TE and re-pellet.
 12. Resuspend the pellet in 50-100 µl TE and plate on a selective plate.
-

Technical Service

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail:
tech_service@invitrogen.com

Japanese Headquarters:

Invitrogen Japan K.K.
Nihonbashi Hama-Cho Park Bldg. 4F
2-35-4, Hama-Cho, Nihonbashi
Tel: 81 3 3663 7972
Fax: 81 3 3663 8242
E-mail: jpinfo@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

MSDS Requests

To request an MSDS, visit our Web site (www.invitrogen.com) and follow the instructions below.

1. On the home page, go to the left-hand column under 'Technical Resources' and select 'MSDS Requests'.
 2. Follow instructions on the page and fill out all the required fields.
 3. To request additional MSDSs, click the 'Add Another' button.
 4. All requests will be faxed unless another method is selected.
 5. When you are finished entering information, click the 'Submit' button. Your MSDS will be sent within 24 hours.
-

continued on next page

Technical Service, continued

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Boder, E. T., and Wittrup, K. D. (1998). Optimal Screening of Surface-Displayed Polypeptide Libraries. *Biotechnology Progress* *14*, 55-62.
- Boder, E. T., and Wittrup, K. D. (1997). Yeast Surface Display for Screening Combinatorial Polypeptide Libraries. *Nature Biotechnology* *15*, 553-557.
- Cappellaro, C., Baldermann, C., Rachel, R., and Tanner, W. (1994). Mating Type-Specific Cell-Cell recognition of *Saccharomyces cerevisiae*: Cell Wall Attachment and Active Sites of α - and α -Agglutinin. *EMBO J.* *13*, 4737-4744.
- Cappellaro, C., Hauser, K., Mersa, V., Watzel, M., Watzel, G., Gruber, C., and Tanner, W. (1991). *Saccharomyces cerevisiae* α - and α -Agglutinin: Characterization of Their Molecular Interaction. *EMBO J.* *10*, 4081-4088.
- Giniger, E., Barnum, S. M., and Ptashne, M. (1985). Specific DNA Binding of GAL4, a Positive Regulatory Protein of Yeast. *Cell* *40*, 767-774.
- Guthrie, C., and Fink, G. R. (1991) Guide to Yeast Genetics and Molecular Biology. In *Methods in Enzymology*, Vol. 194. (J. N. Abelson and M. I. Simon, eds.) Academic Press, San Diego, CA.
- Kieke, M. C., Cho, B. K., Boder, E. T., Kranz, D. M., and Wittrup, K. D. (1997). Isolation of Anti-T Cell Receptor scFv Mutants by Yeast Surface Display. *Protein Engineering* *10*, 1303-1310.
- Lindner, P., Bauer, K., Krebber, A., Nieba, L., Kremmer, E., Krebber, C., Honegger, A., Klinger, B., Mocikat, R., and Pluckthun, A. (1997). Specific Detection of His-tagged Proteins With Recombinant Anti-His Tag scFv-Phosphatase or scFv-Phage Fusions. *BioTechniques* *22*, 140-149.
- Lu, C.-F., Montijn, R. C., Brown, J. L., Klis, F., Kurjan, J., and Bussey, H. (1995). Glycosyl Phosphatidylinositol-Dependent Cross-Linking of Alpha-Agglutinin and Beta 1,6-Glucan in the *Saccharomyces cerevisiae* cell wall. *J. Cell Biol.* *128*, 333-340.
- Roy, A., Lu, C. F., Marykwas, D. L., Lipke, P. N., and Kurjan, J. (1991). The *AGA1* Product Is Involved in Cell Surface Attachment of the *Saccharomyces cerevisiae* Cell Adhesion Glycoprotein α -Agglutinin. *Mol. Cell. Biol.* *11*, 4196-4206.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991). Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. *J. Gen. Virol.* *72*, 1551-1557.
- West, R. W. J., Yocum, R. R., and Ptashne, M. (1984). *Saccharomyces cerevisiae* GAL1-GAL10 Divergent Promoter Region: Location and Function of the Upstream Activator Sequence UAS_G. *Mol. Cell. Biol.* *4*, 2467-2478.