



## Use

When using the pIRESHyg3 Vector, the antibiotic exerts selective pressure on the entire expression cassette; thus, a high dose of antibiotic will select for cells expressing a high level of the gene of interest. This selective pressure also ensures that the expression of the gene of interest will be stable over time in culture. Unless your expression experiments require a pure population of cells, you can use the pool of cells surviving selection instead of isolating and characterizing clonal cell lines. We recommend selecting mammalian cultures in 250–600 µg/ml Hygromycin B (Cat. No. 631309) depending on the cell line. Be sure to establish a kill curve for each cell line and each lot of Hygromycin B to determine optimal selection concentration.

## Location of features

- Human cytomegalovirus (CMV) major immediate early promoter  
Enhancer region: 309–715; TATA Box: 804–810
- T7 RNA polymerase promoter: 863–879
- Multiple cloning site (MCS): 912–1003
- Synthetic intron (IVS): 1002–1297
- Internal ribosome entry site (IRES) from encephalomyocarditis virus (ECMV): 1333–1924
- Hygromycin B phosphotransferase gene: 1948–2979
- SV40 early mRNA polyadenylation signal  
Polyadenylation signals: 3421–3426 & 3450–3455; mRNA 3' ends: 3459 & 3471
- ColE1 origin of replication: 4012–4611
- Ampicillin resistance ( $\beta$ -lactamase) gene:  
Promoter: –35 region: 5703–5698; –10 region: 5680–5675  
 $\beta$ -lactamase coding sequence:  
Start codon: 5633–5631; stop codon: 4775–4773

## Propagation in *E. coli*

- Suitable host strains: DH5 $\alpha$  and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 µg/ml) to *E. coli* hosts.
- *E. coli* replication origin: ColE1
- Copy number: high

## References

1. Jackson, R. J., *et al.* (1990) *Trends Biochem. Sci.* **15**(12):477–483.
2. Jang, S. K., *et al.* (1988) *J. Virol.* **62**(8):2636–2643.
3. Rees, S., *et al.* (1996) *BioTechniques* **20**(1):102–104.
4. Huang, M. T. F. & Gorman, C. M. (1990) *Nucleic Acids Res.* **18**(4):937–947.

**Note:** The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech Laboratories, Inc. This vector has not been completely sequenced.

## Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2008 Clontech Laboratories, Inc.