



pETDuet-1 Vector

TB337 RevA 0903

| | Cat. No. |
|--|-----------------|
| pETDuet-1 DNA | 71146-3 |
| pETDuet-1 sequence landmarks | |
| T7 promoter-1 | 5404–5420 |
| T7 transcription start-1 | 1 |
| His•Tag® coding sequence | 83–100 |
| Multiple cloning sites-1 (<i>Nco</i> I– <i>Afl</i> II) | 69–168 |
| T7 promoter-2 | 214–230 |
| T7 transcription start-2 | 231 |
| Multiple cloning sites-2 (<i>Nde</i> I– <i>Avr</i> II) | 297–438 |
| S•Tag™ coding sequence | 366–410 |
| T7 terminator | 462–509 |
| <i>lacI</i> coding sequence | 3931–5013 |
| pBR322 origin | 2737 |
| <i>bla</i> (<i>Ap</i>) coding sequence | 1119–1976 |
| <i>f1</i> origin | 547–994 |

pETDuet™-1 is designed for the coexpression of two target genes. The vector contains two multiple cloning sites (MCS), each of which is preceded by a T7 promoter/*lac* operator and a ribosome binding site (rbs). The vector also carries the pBR322-derived ColE1 replicon, *lacI* gene and ampicillin resistance gene. This vector can be used in combination with pACYCDuet™-1 (Cat. No. 71147-3) in an appropriate host strain for the coexpression of up to 4 target genes. Genes inserted into MCS1 can be sequenced using the pET Upstream Primer (Cat. No. 69214-3) and DuetDOWN1 Primer (Cat. No. 71179-3). Genes inserted into MCS2 can be sequenced using the DuetUP2 Primer (Cat. No. 71180-3) and T7 Terminator Primer (Cat. No. 69337-3).



pETDuet-1 cloning/expression regions

pETDuet-1 Restriction Sites

| Enzyme | # Sites | Locations | | | |
|----------|---------|----------------|-----------|------|-----------|
| AatII | 1 | 346 | | | |
| AccI | 3 | 138 | 411 | 3026 | |
| AcII | 78 | | | | |
| AfiII | 1 | 163 | | | |
| AfiIII | 2 | 2794 | 4658 | | |
| AhdI | 1 | 1189 | | | |
| AluI | 25 | | | | |
| Alw26I | 7 | 1250 | 2027 | 3153 | 4043 4430 |
| | | 4556 | 4961 | | |
| AlwI | 15 | | | | |
| AlwNI | 1 | 2385 | | | |
| ApaI | 1 | 4455 | | | |
| ApaLI | 4 | 1856 | 2480 | 2980 | 4678 |
| ApoI | 5 | 112 | 384 | 959 | 970 4383 |
| AscI | 1 | 125 | | | |
| AvaI | 1 | 354 | | | |
| AvaII | 5 | 1327 | 1549 | 3510 | 3789 4107 |
| AvrII | 1 | 433 | | | |
| BamHI | 1 | 106 | | | |
| BanI | 9 | 348 | 735 | 1137 | 3889 4019 |
| | | 4738 | 5201 | 5315 | 5336 |
| BanII | 5 | 122 | 705 | 4455 | 5268 5282 |
| BbsI | 3 | 3676 4173 4512 | | | |
| BbvI | 28 | | | | |
| BcgI | 4 | 162 | 1728 | 3223 | 4338 |
| BclI | 1 | 4644 | | | |
| BfaI | 10 | 31 | 415 | 434 | 462 623 |
| | | 1044 | 1379 | 2301 | 3782 3817 |
| BglI | 1 | 1309 | | | |
| BglII | 1 | 305 | | | |
| BpmI | 4 | 1259 | 3273 | 4337 | 4826 |
| Bpu10I | 1 | 3689 | | | |
| Bpu1102I | 1 | 451 | | | |
| BsaAI | 2 | 776 | 3046 | | |
| BsaBI | 3 | 3601 | 5379 | 5389 | |
| BsaHI | 7 | 343 | 1726 | 4020 | 4703 5202 |
| | | 5316 | 5337 | | |
| BsaI | 1 | 1250 | | | |
| BsaJI | 7 | 69 | 433 | 473 | 2634 4023 |
| | | 5215 | 5221 | | |
| BsaWI | 7 | 528 | 1480 | 2441 | 2588 3605 |
| | | 3836 | 4339 | | |
| BsgI | 3 | 3640 4613 4813 | | | |
| BsiEI | 9 | 153 | 199 | 325 | 337 1559 |
| | | 1708 | 2460 | 2884 | 3879 |
| BsiHKAI | 8 | 122 | 1775 | 1860 | 2484 2984 |
| | | 3808 | 4682 5166 | | |
| BsII | 22 | | | | |
| BsmBI | 2 | 3153 | 4043 | | |
| BsmFI | 3 | 557 | 3523 | 5197 | |
| Bsp1286I | 12 | | | | |
| BspEI | 2 | 528 | 3605 | | |
| BspLU11I | 1 | 2794 | | | |
| BspMI | 1 | 124 | | | |
| BsrBI | 5 | 13 | 243 | 632 | 2031 2865 |
| BsrDI | 4 | 1250 | 1424 | 4251 | 4617 |
| BsrFI | 6 | 324 | 671 | 1269 | 4972 5339 |
| | | 5348 | | | |
| BsrGI | 1 | 190 | | | |
| BsrI | 24 | | | | |
| BssHII | 2 | 125 | 4247 | | |
| BssSI | 2 | 1853 | 2621 | | |
| Bst1107I | 1 | 3027 | | | |
| BstEII | 1 | 4476 | | | |
| BstXI | 3 | 4612 | 4735 | 4864 | |
| BstYI | 11 | | | | |
| Cac8I | 36 | | | | |
| Clal | 1 | 5383 | | | |
| CviJI | 84 | | | | |
| Ddel | 10 | 262 | 451 | 1146 | 1686 2111 |
| | | 2520 | 2987 | 3527 | 3689 4087 |
| DpnI | 28 | | | | |

| Enzyme | # Sites | Locations | | | | |
|----------|---------|-------------|------|------|-----------|--|
| DraI | 3 | 1055 | 1074 | 1766 | | |
| DraIII | 1 | 779 | | | | |
| DrdI | 3 | 823 | 2692 | 3107 | | |
| DsaI | 2 | 69 | 5221 | | | |
| EaeI | 8 | 150 | 196 | 322 | 326 1577 | |
| | | 3984 | 5218 | 5350 | | |
| EagI | 3 | 150 196 322 | | | | |
| EarI | 3 | 1984 | 2911 | 5041 | | |
| Ecl136II | 1 | 120 | | | | |
| Eco47III | 2 | 3544 | 5257 | | | |
| Eco57I | 2 | 1856 2252 | | | | |
| EcoNI | 1 | 5126 | | | | |
| EcoO109I | 3 | 478 | 3789 | 5226 | | |
| EcoRI | 1 | 112 | | | | |
| EcoRII | 8 | 102 | 2633 | 2646 | 2767 4022 | |
| | | 4079 | 4619 | 4934 | | |
| EcoRV | 1 | 319 | | | | |
| EheI | 4 | 4021 | 5203 | 5317 | 5338 | |
| FauI | 16 | | | | | |
| Fnu4HI | 48 | | | | | |
| FokI | 10 | 1155 | 1336 | 1623 | 3108 3249 | |
| | | 3435 | 3513 | 3575 | 4603 4612 | |
| FseI | 1 | 328 | | | | |
| FspI | 1 | 1411 | | | | |
| HaeII | 13 | | | | | |
| HaeIII | 24 | | | | | |
| HgaI | 12 | | | | | |
| HhaI | 44 | | | | | |
| HincII | 2 | 139 | 4156 | | | |
| HindIII | 1 | 143 | | | | |
| Hinfl | 16 | | | | | |
| HpaI | 1 | 4156 | | | | |
| HphI | 18 | | | | | |
| KpnI | 1 | 352 | | | | |
| MaeIII | 17 | | | | | |
| MboII | 13 | | | | | |
| MluI | 1 | 4658 | | | | |
| MnlI | 25 | | | | | |
| MseI | 33 | | | | | |
| MsiI | 9 | 1441 | 1600 | 1959 | 3225 3616 | |
| | | 3811 | 4292 | 4322 | 4610 | |
| MspA1I | 10 | 375 | 450 | 1824 | 2211 2456 | |
| | | 3088 | 3207 | 3969 | 4062 4632 | |
| MspI | 29 | | | | | |
| MunI | 1 | 311 | | | | |
| MwoI | 34 | | | | | |
| NarI | 4 | 4020 | 5202 | 5316 | 5337 | |
| NciI | 12 | | | | | |
| NcoI | 1 | 69 | | | | |
| NdeI | 1 | 298 | | | | |
| NgoAIV | 3 | 324 | 671 | 5348 | | |
| NlaIII | 26 | | | | | |
| NlaIV | 22 | | | | | |
| NotI | 1 | 150 | | | | |
| NspI | 4 | 2798 | 3165 | 3457 | 5191 | |
| Pacl | 1 | 429 | | | | |
| PfIMI | 2 | 401 5083 | | | | |
| PleI | 12 | | | | | |
| Psp1406I | 5 | 989 | 1415 | 1788 | 3471 4998 | |
| Psp5II | 1 | 3789 | | | | |
| PstI | 1 | 135 | | | | |
| PvuI | 2 | 337 | 1559 | | | |
| PvuII | 3 | 3207 | 3969 | 4062 | | |
| RcaI | 4 | 1993 | 2025 | 2074 | 5260 | |
| RsaI | 5 | 192 | 350 | 1669 | 2992 4515 | |
| SacI | 1 | 122 | | | | |
| Sall | 1 | 137 | | | | |
| SapI | 1 | 2911 | | | | |
| Sau3AI | 28 | | | | | |
| Sau96I | 15 | | | | | |
| Scal | 1 | 1669 | | | | |
| ScrFI | 20 | | | | | |

| Enzyme | # Sites | Locations | | | |
|----------|---------|-----------|------|------|-----------|
| SfaNI | 20 | | | | |
| SfcI | 7 | 131 | 226 | 553 | 1430 2338 |
| | | 2529 | 5416 | | |
| SgfI | 1 | 337 | | | |
| SgrAI | 1 | 5339 | | | |
| SphI | 1 | 5191 | | | |
| Sse8387I | 1 | 135 | | | |
| SspI | 1 | 984 | | | |
| StyI | 3 | 69 | 433 | 473 | |
| TaiI | 15 | | | | |
| TaqI | 16 | | | | |
| TfiI | 4 | 2820 | 3241 | 3745 | 3980 |
| Thal | 34 | | | | |
| TseI | 28 | | | | |
| Tsp45I | 7 | 598 | 1445 | 1656 | 3046 3141 |
| | | 3354 | 4476 | | |
| Tsp509I | 22 | | | | |
| TspRI | 13 | | | | |
| Tth111I | 1 | 3052 | | | |
| VspI | 5 | 213 | 1361 | 3916 | 3975 5403 |
| XbaI | 1 | 30 | | | |
| XcmI | 3 | 4273 | 4291 | 4807 | |
| XhoI | 1 | 354 | | | |
| XmnI | 2 | 1788 | 3240 | | |

Enzymes that do not cut pETDuet-1:

| | | | | | |
|-------|-------|--------|-------|------|-------|
| BseRI | BsmI | Bsu36I | MscI | NheI | NruI |
| NsiI | NspV | PinAI | PmeI | PmlI | PshAI |
| RsrII | SacII | SanDI | SexAI | SfiI | Smal |
| SnaBI | SpeI | SrfI | StuI | SunI | Swal |