

# pRL-TK Vector



Technical Bulletin No. 240

INSTRUCTIONS FOR USE OF PRODUCT E2241. PLEASE DISCARD PREVIOUS VERSIONS.

All technical literature is available on the Internet at [www.promega.com](http://www.promega.com)

Please visit the web site to verify that you are using the most current version of this Technical Bulletin.

I. Description .....	1
II. Product Components .....	2
III. Features of the pRL-TK Vector .....	3
A. TK Promoter Region .....	3
B. Chimeric Intron .....	3
C. T7 Promoter .....	3
D. <i>Renilla</i> Luciferase Reporter Gene ( <i>Rluc</i> ) .....	3
E. SV40 Late Polyadenylation Signal .....	3
IV. Transfection of Mammalian Cells with pRL-TK .....	4
V. pRL-TK Vector Restriction Sites and Sequence .....	4
A. pRL-TK Vector Restriction Sites .....	4
B. pRL-TK Vector Sequence .....	6
VI. Related Products .....	9
VII. References .....	9

## I. Description

The pRL-TK Vector<sup>(a,b)</sup> (Figure 1) is intended for use as an internal control reporter and may be used in combination with any experimental reporter vector to cotransfect mammalian cells. All of Promega's pRL Reporter Vectors contain a cDNA (*Rluc*) encoding *Renilla* luciferase, which was originally cloned from the marine organism *Renilla reniformis* (sea pansy; 1). As described below, the *Renilla* luciferase cDNA<sup>(b)</sup> contained within the pRL Vectors has been modified slightly to provide greater utility.

The pRL-TK Vector contains the herpes simplex virus thymidine kinase (HSV-TK) promoter to provide low to moderate levels of *Renilla* luciferase expression in co-transfected mammalian cells. *Renilla* luciferase is a 36kDa monomeric protein that does not require post-translational modification for activity (2). Therefore, like firefly luciferase, the enzyme may function as a genetic reporter immediately following translation. For information about the use of this plasmid in conjunction with a reporter vector containing the firefly luciferase gene, refer to the *Dual-Luciferase<sup>®</sup> Reporter Assay System<sup>(c,d)</sup> Technical Manual* (#TM040).

The pRL Vectors are isolated from a *dam*<sup>-</sup>/*dcm*<sup>-</sup> *E. coli* K host strain, allowing digestion with restriction enzymes that are sensitive to *dam* and *dcm* methylation.

The GenBank<sup>®</sup>/EMBL Accession number for the pRL-TK Vector is AF025846.




## II. Product Components

Product	Size	Cat.#
pRL-TK Vector	20µg	E2241

All pRL Vectors are supplied in TE buffer (pH 7.4) and are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain vector and are not competent cells.

**Storage Conditions:** Store vector DNA at  $-20^{\circ}\text{C}$  and the glycerol stock of JM109 cells at  $-70^{\circ}\text{C}$ .




**Figure 1. The pRL-TK Vector circle map and sequence reference points.**

### Sequence reference points:

HSV-TK promoter	7-759
Chimeric intron	826-962
T7 RNA polymerase promoter (-17 to +2)	1006-1024
T7 RNA polymerase transcription initiation site	1023
<i>RLuc</i> reporter gene	1034-1969
SV40 late polyadenylation signal	2011-2212
$\beta$ -lactamase ( <i>Amp<sup>r</sup></i> ) coding region	2359-3219

### In addition:

-  indicates the position of the intron.
- *RLuc* is the cDNA encoding the *Renilla* luciferase enzyme.
- *Amp<sup>r</sup>* indicates the gene encoding ampicillin resistance in *E. coli*.
- ori is the origin of replication in *E. coli*.
- The arrows within the *RLuc* and *Amp<sup>r</sup>* genes indicate the direction of transcription. Restriction sites shown in parentheses are not unique sites.

### III. Features of the pRL-TK Vector

#### A. TK Promoter Region

The pRL-TK Vector contains the herpes simplex virus thymidine kinase promoter region upstream of *Rluc*. The HSV-TK promoter provides low-level, constitutive expression in cells of both embryonal and mature mammalian tissues (3,4).

#### B. Chimeric Intron

Downstream of the TK promoter region of the pRL-TK Vector is a chimeric intron comprised of the 5'-donor splice site from the first intron of the human  $\beta$ -globin gene, and the branch and 3'-acceptor splice site from an intron preceding an immunoglobulin gene heavy chain variable region (5). The sequences of the donor and acceptor splice sites, along with the branchpoint site, have been modified to match the consensus sequences for optimal splicing (6).

Transfection studies have demonstrated that the presence of an intron flanking a cDNA insert frequently increases the level of gene expression (7–10). In the pRL-SV40 Vector, the intron is positioned 5' to *Rluc* to minimize the utilization of cryptic 5'-donor splice sites that may reside within the reporter gene sequence (11).

#### C. T7 Promoter

A T7 promoter is located downstream of the chimeric intron, immediately preceding the *Rluc* reporter gene. This T7 promoter can be used to synthesize RNA transcripts in vitro using T7 RNA Polymerase (Cat.# P2075). T7 RNA Polymerase can also be used to synthesize active *Renilla* luciferase in a cell-free coupled eukaryotic in vitro transcription/translation reaction (e.g., Promega's TNT® Reticulocyte Lysate [Cat.# L4610]<sup>(c,e,f,g)</sup>, Wheat Germ Extract<sup>(c,e,f,g)</sup> [Cat.# L4140], or TNT® T7 Quick Coupled Transcription/Translation<sup>(c,e,f,g,h)</sup> [Cat.# L1170] Systems).

#### D. *Renilla* Luciferase Reporter Gene (*Rluc*)

The *Renilla* luciferase cDNA inserted into all of the pRL Vectors is derived from the anthozoan coelenterate *Renilla reniformis* (1) but contains nucleotide changes that were engineered during the construction of the individual vectors. The following bases were altered in the pRL-TK Vector: base 1264 (T→C), to eliminate an internal *Bgl* II site; base 1807 (T→C), to eliminate internal *Bam*H I site; base 1840 (C→T), to eliminate internal *Nar* I, *Kas* K, *Ban* K and *Acy* I sites. These nucleotide substitutions do not alter the amino acid sequence of the encoded *Renilla* luciferase reporter enzyme.

#### E. SV40 Late Polyadenylation Signal

Polyadenylation signals cause the termination of transcription by RNA polymerase II and signal the addition of approximately 200–250 adenosine residues to the 3'-end of the RNA transcript (12). Polyadenylation has been shown to enhance RNA stability and translation (13,14). The late SV40 polyadenylation signal, which is extremely efficient and has been shown to increase the steady-state level of RNA approximately five-fold more than the early SV40 polyadenylation signal (15), has been positioned 3' to the *Rluc* gene in the pRL-TK Vector to increase the level of *Renilla* luciferase expression.



**Due to** sequence differences between the T7 Promoter Primer offered by Promega (Cat.# Q5021) and the T7 promoter used in the pRL family of *Renilla* luciferase co-reporter vectors, Cat.# Q5021 cannot be used for sequencing with this vector.

#### IV. Transfection of Mammalian Cells with pRL-TK

The pRL-TK Vector may be used in combination with any experimental reporter vector to cotransfect mammalian cells. However, it is important to realize that *trans* effects between promoters on cotransfected plasmids can potentially affect reporter gene expression (16). This is primarily of concern when either the control or experimental reporter vector, or both, contain very strong promoter/enhancer elements. The occurrence and magnitude of such effects will depend on several factors: i) the combination and activities of the genetic regulatory elements present on the cotransfected vectors, ii) the relative ratio of experimental vector to control vector introduced into the cells, and iii) the type of cell transfected.

To help ensure independent genetic expression between experimental and control reporter genes, preliminary cotransfection experiments should be performed to optimize both the amount of vector DNA and the ratio of the coreporter vectors added to the transfection mixture. Similar to the firefly luciferase assay, the *Renilla* luciferase assay is extremely sensitive, providing accurate measurement of  $\leq 10$  femtograms of *Renilla* luciferase, with linearity over 7 orders of enzyme concentration. Therefore, it is possible to use relatively small quantities of pRL-TK Vector to provide low-level, constitutive coexpression of *Renilla* luciferase control activity. Ratios of 10:1 (or greater) for experimental vector:pRL-TK Vector combinations may aid greatly in suppressing the occurrence of *trans* effects between promoter elements.

The pRL-TK Vector can be used for both transient and stable expression of genes. For stable expression, the pRL-TK Vector must be cotransfected with an expression vector containing a selectable gene in mammalian cells. Transfection of DNA into mammalian cells may be mediated by cationic lipids (17,18), calcium phosphate (19,20), DEAE-Dextran (21–23), polybrene-DMSO (24,25) or electroporation (26,27).

Transfection systems based on cationic lipid compounds (TransFast™ Reagent<sup>(i)</sup>, Tfx™ Reagents<sup>(i)</sup> and Transfectam® Reagent<sup>(k)</sup>), calcium phosphate and DEAE-Dextran are available from Promega. For more information and a protocol for the Transfectam® Reagent, please request the *Transfectam® Reagent Technical Bulletin* (#TB116) and for TransFast™ Reagent, please request the *TransFast™ Transfection Reagent Technical Bulletin* (#TB260). Protocols for the use of the Tfx™ Reagents can be found in the *Tfx™-10*, *Tfx™-20* and *Tfx™-50 Reagent Technical Bulletin* (#TB216). For transfection procedures using calcium phosphate or DEAE-Dextran, please request the *ProFection® Mammalian Transfection Systems Technical Manual* (#TM012).

#### V. pRL-TK Vector Restriction Sites and Sequence

##### A. pRL-TK Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch Office or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526.

**Table 1. Restriction Enzymes That Cut the pRL-TK Vector Between 1 and 5 Times.**

<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>	<b>Enzyme</b>	<b># of Sites</b>	<b>Location</b>
<b>Acc I</b>	1	342	<i>Drd I</i>	2	781, 3932
<b>Acy I</b>	2	290, 2606	<i>Dsa I</i>	5	129, 180, 296, 330 543
<i>Afl II</i>	4	36, 160, 792, 989	<i>Eae I</i>	4	126, 1384, 1978, 2753
<i>Afl III</i>	2	713, 1216	<i>Eag I</i>	1	1978
<b>Alw44 I</b>	2	2474, 3720	<i>Ear I</i>	3	137, 1204, 2347
<i>AlwN I</i>	2	357, 3625	<b>EcIHK I</b>	1	3146
<i>AspH I</i>	3	2478, 2563, 3724	<b>Eco52 I</b>	1	1978
<b>Ava I</b>	3	108, 229, 282	<b>EcoR I</b>	1	649
<i>Avr II</i>	1	322	<i>Ehe I</i>	1	291
<b>Bal I</b>	1	128	<i>Fsp I</i>	1	2923
<b>BamH I</b>	1	2223	<b>Hae II</b>	3	150, 293, 3794
<b>Ban I</b>	4	289, 915, 1838, 3193	<b>Hinc II</b>	1	2121
<b>Ban II</b>	1	313	<i>Hind II</i>	1	2121
<i>Bbe I</i>	1	293	<b>Hind III</b>	1	760
<i>Bbs I</i>	3	6, 900, 1874	<b>Hpa I</b>	1	2121
<b>Bcl I</b>	2	1318, 1527	<b>Hsp92 I</b>	2	290, 2606
<b>Bgl I</b>	3	298, 301, 3028	<i>Kas I</i>	1	289
<b>Bgl II</b>	1	1	<i>Mlu I</i>	1	713
<i>Bsa I</i>	3	386, 854, 3080	<b>Nar I</b>	1	290
<b>BsaO I</b>	4	1981, 2628, 2777, 3700	<b>Nco I</b>	1	330
<i>BsaA I</i>	2	346, 1766	<b>Nhe I</b>	1	1024
<i>BsaB I</i>	1	2222	<b>Not I</b>	1	1978
<i>BsaH I</i>	2	290, 2606	<i>Nsp I</i>	3	599, 1160, 1220
<b>BsaM I</b>	2	2042, 2135	<i>PaeR7 I</i>	2	882, 4040
<i>Bsm I</i>	2	2042, 2135	<i>PpuM I</i>	1	319
<i>BspH I</i>	3	1602, 2306, 3314	<i>Psp5 II</i>	1	319
<i>BspM I</i>	1	816	<i>PspA I</i>	1	282
<b>BsrBR I</b>	1	2222	<b>Pst I</b>	3	509, 746, 802
<i>BsrG I</i>	1	1732	<b>Pvu I</b>	1	2777
<b>BssH II</b>	1	241	<b>Pvu II</b>	1	531
<i>BssS I</i>	3	1692, 2477, 3861	<b>Rsa I</b>	4	473, 1002, 1734, 2665
<i>Bst1107 I</i>	1	343	<b>Sac II</b>	1	299
<b>Bst98 I</b>	4	36, 160, 792, 989	<b>Sca I</b>	2	1002, 2665
<b>BstZ I</b>	1	1978	<b>Sma I</b>	1	284
<i>Cfr10 I</i>	1	3061	<b>Ssp I</b>	1	2341
<b>Cla I</b>	1	2216	<b>Sty I</b>	2	322, 330
<b>Csp45 I</b>	2	653, 1040	<b>Vsp I</b>	2	1134, 2971
<b>Dde I</b>	4	2645, 3185, 3351, 3760	<b>Xba I</b>	1	1971
<b>Dra I</b>	4	2182, 2568, 3260, 3279	<i>Xcm I</i>	1	1683
<i>Dra II</i>	1	319	<b>Xma I</b>	1	282
			<b>Xmn I</b>	2	1568, 2546

**Note:** The enzymes listed in boldface type are available from Promega.

**Table 2. Restriction Enzymes That Do Not Cut the pRL-TK Vector.**

<i>Aat</i> II	<i>Bpu</i> 1102 I	<b><i>Eco</i>ICR I</b>	<i>Nsi</i> I	<b><i>Sal</i> I</b>	<i>Sse</i> 8387 I
<b><i>Acc</i>B7 I</b>	<i>Bsp</i> 120 I	<i>Eco</i> N I	<i>Pac</i> I	<b><i>Sfi</i> I</b>	<b><i>Stu</i> I</b>
<b><i>Acc</i> III</b>	<b><i>Bst</i>E II</b>	<b><i>Eco</i>R V</b>	<i>Pfl</i> M I	<b><i>Sgf</i>I<sup>(1)</sup></b>	<i>Swa</i> I
<b><i>Acc</i>65I</b>	<b><i>Bst</i>X I</b>	<i>Fse</i> I	<i>Pin</i> A I	<i>Sgr</i> A I	<b><i>Tth</i>111 I</b>
<b><i>Age</i> I</b>	<b><i>Bsu</i>36 I</b>	<b><i>I-Ppo</i> I</b>	<i>Pme</i> I	<b><i>Sna</i>B I</b>	<b><i>Xho</i> I</b>
<b><i>Apa</i> I</b>	<b><i>Csp</i> I</b>	<b><i>Kpn</i> I</b>	<i>Pml</i> I	<b><i>Spe</i> I</b>	
<i>Asc</i> I	<i>Dra</i> III	<b><i>Nae</i> I</b>	<i>Ppu</i> 10 I	<b><i>Sph</i> I</b>	
<i>Bbr</i> P I	<b><i>Eco</i>47 III</b>	<b><i>Nde</i> I</b>	<i>Psh</i> A I	<i>Spl</i> I	
<b><i>Bbu</i> I</b>	<i>Eco</i> 72 I	<b><i>Ngo</i>MIV</b>	<i>Rsr</i> II	<i>Srf</i> I	
<i>Blp</i> I	<i>Eco</i> 81 I	<b><i>Nru</i> I</b>	<b><i>Sac</i> I</b>		

**Note:** The enzymes listed in boldface type are available from Promega.

**Table 3. Restriction Enzymes That Cut the pRL-TK Vector 6 or More Times.**

<i>Aci</i> I	<b><i>Bsr</i>S I</b>	<b><i>Fok</i> I</b>	<i>Mae</i> I	<b><i>Msp</i>A1 I</b>	<i>Scr</i> F I
<b><i>Alu</i> I</b>	<b><i>Bst</i>71 I</b>	<b><i>Hae</i> III</b>	<i>Mae</i> II	<b><i>Nci</i> I</b>	<i>Sfa</i> N I
<b><i>Alw</i>26 I</b>	<b><i>Bst</i>O I</b>	<i>Hga</i> I	<i>Mae</i> III	<b><i>Nde</i> II</b>	<b><i>Sin</i> I</b>
<b><i>Ava</i> II</b>	<i>Bst</i> U I	<b><i>Hha</i> I</b>	<b><i>Mbo</i> I</b>	<i>Nla</i> III	<b><i>Taq</i> I</b>
<i>Bbv</i> I	<b><i>Cfo</i> I</b>	<b><i>Hinf</i> I</b>	<b><i>Mbo</i> II</b>	<i>Nla</i> IV	<i>Tfi</i> I
<i>Bsa</i> J I	<b><i>Dpn</i> I</b>	<b><i>Hpa</i> II</b>	<i>Mnl</i> I	<i>Ple</i> I	<b><i>Tru</i>9 I</b>
<b><i>Bsp</i>1286 I</b>	<i>Dpn</i> II	<i>Hph</i> I	<i>Mse</i> I 20	<b><i>Sau</i>3A I</b>	<b><i>Xho</i> II</b>
<i>Bsr</i> I	<i>Fnu</i> 4H I	<b><i>Hsp</i>92 II</b>	<b><i>Msp</i> I</b>	<b><i>Sau</i>96 I</b>	

**Note:** The enzymes listed in boldface type are available from Promega.

## B. pRL-TK Vector Sequence

The sequence shown corresponds to the mRNA synthesized from the *Renilla* luciferase gene from the TK promoter. The vector sequence is also available on the Internet at [www.promega.com/vectors](http://www.promega.com/vectors). The GenBank®/EMBL Accession Number for the pRL-TK Vector is AF025846.

```

1  AGATCTAAAT GAGTCTTCGG ACCTCGCGGG GCGCGCTTAA GCGGTGGTTA
51  GGGTTTGTCT GACGCGGGGG GAGGGGGAAG GAACGAAACA CTCTCATTCG
101 GAGGCGGCTC GGGGTTTGGT CTTGGTGGCC ACGGGCACGC AGAAGAGCGC
151 CGCGATCCTC TTAAGCACCC CCCC GCCCTC CGTGGAGGCG GGGGTTTGGT
201 CGGCGGGTGG TAACTGGCGG GCCGCTGACT CGGGCGGGTC GCGCGCCCCA
251 GAGTGTGACC TTTTCGGTCT GCTCGCAGAC CCCC GGGCGG CGCCGCCGCG
301 GCGGCGACGG GCTCGCTGGG TCCTAGGCTC CATGGGGACC GTATACGTGG
351 ACAGGCTCTG GAGCATCCGC ACGACTGCGG TGATATTACC GGAGACCTTC
401 TGCGGGACGA GCCGGGTCAC GCGGCTGACG CGGAGCGTCC GTTGGGCGAC
451 AAACACCAGG ACGGGGCACA GGTACACTAT CTTGTCAACC GGAGGCGCGA
501 GGGACTGCAG GAGCTTCAGG GAGTGGCGCA GCTGCTTCAT CCCC GTGGCC
551 CGTTGCTCGC GTTTGCTGGC GGTGTCCCCG GAAGAAATAT ATTTGCATGT
601 CTTTAGTTCT ATGATGACAC AAACCCCGCC CAGCGTCTTG TCATTGGCGA
651 ATTCGAACAC GCAGATGCAG TCGGGGCGGC GCGGTCCCAG GTCCACTTCG

```

**B. pRL-TK Vector Sequence (continued)**

```

701 CATATTAAGG TGACGCGTGT GGCCTCGAAC ACCGAGCGAC CCTGCAGCGA
751 CCCGCTTAAA AGCTTGATTG TTCTGACACA ACAGTCTCGA ACTTAAGCTG
801 CAGAAGTTGG TCGTGAGGCA CTGGGCAGGT AAGTATCAAG GTTACAAGAC
851 AGGTTTAAGG AGACCAATAG AAAGTGGGCT TGTCGAGACA GAGAAGACTC
901 TTGCGTTTCT GATAGGCACC TATTGGTCTT ACTGACATCC ACTTTGCCTT
951 TCTCTCCACA GGTGTCCACT CCCAGTTCOA TTACAGCTCT TAAGGCTAGA
1001 GTACTTAATA CGACTCACTA TAGGCTAGCC ACCATGACTT CGAAAGTTTA
1051 TGATCCAGAA CAAAGGAAAC GGATGATAAC TGGTCCGCAG TGGTGGGCCA
1101 GATGTAAACA AATGAATGTT CTTGATTCAT TTATTAATTA TTATGATTCA
1151 GAAAAACATG CAGAAAATGC TGTATTTTTT TTACATGGTA ACGCGGCCTC
1201 TTCTTATTTA TGGCGACATG TTGTGCCACA TATTGAGCCA GTAGCGCGGT
1251 GTATTATACC AGACCTTATT GGTATGGGCA AATCAGGCAA ATCTGGTAAT
1301 GGTTCTTATA GGTTACTTGA TCATTACAAA TATCTTACTG CATGGTTTGA
1351 ACTTCTTAAT TTACCAAAGA AGATCATTTT TGTCGGCCAT GATTGGGGTG
1401 CTTGTTTGGC ATTTCAATTAT AGCTATGAGC ATCAAGATAA GATCAAAGCA
1451 ATAGTTCACG CTGAAAGTGT AGTAGATGTG ATTGAATCAT GGGATGAATG
1501 GCCTGATATT GAAGAAGATA TTGCGTTGAT CAAATCTGAA GAAGGAGAAA
1551 AAATGGTTTT GGAGAATAAC TTCTTCGTGG AAACCATGTT GCCATCAAAA
1601 ATCATGAGAA AGTTAGAACC AGAAGAATTT GCAGCATATC TTGAACCATT
1651 CAAAGAGAAA GGTGAAGTTC GTCGTCCAAC ATTATCATGG CCTCGTGAAA
1701 TCCCGTTAGT AAAAGGTGGT AAACCTGACG TTGTACAAAT TGTTAGGAAT
1751 TATAATGCTT ATCTACGTGC AAGTGATGAT TTACCAAAAA TGTTTATTGA
1801 ATCGGACCCA GGATTCTTTT CCAATGCTAT TGTTGAAGGT GCCAAGAAGT
1851 TTCCTAATAC TGAATTTGTC AAAGTAAAAG GTCTTCATTT TTCGCAAGAA
1901 GATGCACCTG ATGAAATGGG AAAATATATC AAATCGTTCG TTGAGCGAGT
1951 TCTCAAAAAT GAACAATAAT TCTAGAGCGG CCGCTTCGAG CAGACATGAT
2001 AAGATACATT GATGAGTTTG GACAAACCAC AACTAGAATG CAGTGAAAAA
2051 AATGCTTTAT TTGTGAAATT TGTGATGCTA TTGCTTTATT TGTAACCATT
2101 ATAAGCTGCA ATAAACAAGT TAACAACAAC AATTGCATTC ATTTTATGTT
2151 TCAGGTTTCA GGGGAGGTGT GGGAGGTTTT TTAAAGCAAG TAAAACCTCT
2201 ACAAATGTGG TAAAATCGAT AAGGATCCAG GTGGCACTTT TCGGGGAAAT
2251 GTGCGCGGAA CCCCTATTTG TTTATTTTTC TAAATACATT CAAATATGTA
2301 TCCGCTCATG AGACAATAAC CCTGATAAAT GCTTCAATAA TATTGAAAAA
2351 GGAAGAGTAT GAGTATTCAA CATTTCCGTG TCGCCCTTAT TCCCTTTTTT
2401 GCGGCATTTT GCCTTCCTGT TTTTGCTCAC CCAGAAACGC TGGTGAAAGT
2451 AAAAGATGCT GAAGATCAGT TGGGTGCACG AGTGGGTTAC ATCGAACTGG
2501 ATCTCAACAG CGGTAAGATC CTTGAGAGTT TTCGCCCGA AGAACGTTTT

```



B. pRL-TK Vector Sequence (continued)

```
2551 CCAATGATGA GCACTTTTAA AGTTCTGCTA TGTGGCGCGG TATTATCCCG
2601 TATTGACGCC GGGCAAGAGC AACTCGGTCG CCGCATAAC TATTCTCAGA
2651 ATGACTTGGT TGAGTACTCA CCAGTCACAG AAAAGCATCT TACGGATGGC
2701 ATGACAGTAA GAGAATTATG CAGTGCTGCC ATAACCATGA GTGATAACAC
2751 TGC GGCCAAC TTACTTCTGA CAACGATCGG AGGACCGAAG GAGCTAACCG
2801 CTTTTTTGCA CAACATGGGG GATCATGTAA CTCGCCTTGA TCGTTGGGAA
2851 CCGGAGCTGA ATGAAGCCAT ACCAAACGAC GAGCGTGACA CCACGATGCC
2901 TG TAGCAATG GCAACAACGT TGCGCAAAC ATTAAC TGGC GAACTACTTA
2951 CTCTAGCTTC CCGGCAACAA TTAATAGACT GGATGGAGGC GGATAAAGTT
3001 GCAGGACCAC TTCTGCGCTC GGCCCTTCCG GCTGGCTGGT TTATTGCTGA
3051 TAAATCTGGA GCCGGTGAGC GTGGGTCTCG CGGTATCATT GCAGCACTGG
3101 GGCCAGATGG TAAGCCCTCC CGTATCGTAG TTATCTACAC GACGGGGAGT
3151 CAGGCAACTA TGGATGAACG AAATAGACAG ATCGCTGAGA TAGGTGCCTC
3201 ACTGATTAAG CATTGGTAAC TGTCAGACCA AGTTTACTCA TATATACTTT
3251 AGATTGATTT AAAACTTCAT TTTTAATTTA AAAGGATCTA GGTGAAGATC
3301 CTTTTTGATA ATCTCATGAC CAAAATCCCT TAACGTGAGT TTTCTTCCA
3351 CTGAGCGTCA GACCCCGTAG AAAAGATCAA AGGATCTTCT TGAGATCCTT
3401 TTTTTCTGCG CGTAATCTGC TGCTTGCAA CAAAAAACC ACCGCTACCA
3451 GCGGTGGTTT GTTTGCCGGA TCAAGAGCTA CCAACTCTTT TTCCGAAGGT
3501 AACTGGCTTC AGCAGAGCGC AGATACCAA TACTGTTCTT CTAGTGTAGC
3551 CGTAGTTAGG CCACCACTTC AAGAACTCTG TAGCACCGCC TACATACCTC
3601 GCTCTGCTAA TCCTGTTACC AGTGGCTGCT GCCAGTGGCG ATAAGTCGTG
3651 TCTTACCGGG TTGGA CTCAA GACGATAGTT ACCGGATAAG GCGCAGCGGT
3701 CGGGCTGAAC GGGGGGTTTCG TGCACACAGC CCAGCTTGGG GCGAACGACC
3751 TACACCGAAC TGAGATACCT ACAGCGTGAG CTATGAGAAA GCGCCACGCT
3801 TCCCGAAGGG AGAAAGGCGG ACAGGTATCC GGTAAGCGGC AGGGTCGGAA
3851 CAGGAGAGCG CACGAGGGAG CTTCCAGGGG GAAACGCCTG GTATCTTTAT
3901 AGTCCTGTCG GGTTTCGCCA CCTCTGACTT GAGCGTCGAT TTTTGTGATG
3951 CTCGTCAGGG GGGCGGAGCC TATGGAAAAA CGCCAGCAAC GCGGCCTTTT
4001 TACGGTTCCT GGCCTTTTGC TGGCCTTTTG CTCACATGGC TCGAC
```



## VI. Related Products

### pRL Family of *Renilla* Luciferase Vectors for Co-Reporter Applications

Product	Size	Cat.#
pRL-CMV Vector <sup>(a,b,m)</sup>	20µg	E2261
pRL-SV40 Vector <sup>(a,b)</sup>	20µg	E2231
pRL-null Vector <sup>(a,b)</sup>	20µg	E2271

To inquire about the availability of bulk packaging and pricing for pRL Vectors, please contact Promega. For inquiries on the availability of new promoter variations within the pRL family of co-reporter vectors, contact Technical Services or visit our web site at [www.promega.com](http://www.promega.com).

### pGL3 Vectors

Product	Size	Cat.#
pGL3-Control Vector <sup>(a,f,n)</sup>	20µg	E1741
pGL3-Basic Vector <sup>(a,f,n)</sup>	20µg	E1751
pGL3-Promoter Vector <sup>(a,f,n)</sup>	20µg	E1761
pGL3-Enhancer Vector <sup>(a,f,n)</sup>	20µg	E1771

### Luciferase Assay Systems

Product	Size	Cat.#
Dual-Luciferase <sup>®</sup> Reporter Assay System	100 assays	E1910
Dual-Luciferase <sup>®</sup> Reporter Assay System 10-pack <sup>(c,d)</sup>	1,000 assays	E1960
Dual-Luciferase <sup>®</sup> Reporter 1000 Assay System <sup>(c,d)</sup>	1,000 assays	E1980

### Transfection Systems

Product	Size	Cat.#
TransFast <sup>™</sup> Transfection Reagents	1.2mg	E2431
Transfectam <sup>®</sup> Reagent for the Transfection of Eukaryotic Cells	1mg	E1231
	0.5mg	E1232
Tfx <sup>™</sup> -50 Reagent	2.1mg	E1811
Tfx <sup>™</sup> -20 Reagent	4.8mg	E2391
Tfx <sup>™</sup> -10 Reagent	9.3mg	E2381
Tfx <sup>™</sup> Reagents Transfection Trio	5.4mg	E2400
ProFection <sup>®</sup> Mammalian Transfection System —Calcium Phosphate	80 transfections	E1200
ProFection <sup>®</sup> Mammalian Transfection System —DEAE-Dextran	80 transfections	E1210

## VII. References

- Lorenz, W.W. *et al.* (1991) Isolation and expression of a cDNA encoding *Renilla reniformis* luciferase. *Proc. Natl. Acad. Sci. USA* **88**, 4438–42.
- Matthews, J.C. *et al.* (1977) Substrate and substrate analogue binding properties of *Renilla* luciferase. *Biochemistry* **16**, 85–91.
- Wagner, E.F. *et al.* (1985) Transfer of genes into embryonal carcinoma cells by retrovirus infection: efficient expression from an internal promoter. *EMBO J.* **4**, 663–6.
- Stewart, C.L. *et al.* (1987) Expression of retroviral vectors in transgenic mice obtained by embryo infection. *EMBO J.* **6**, 383–8.
- Bothwell, A.L.M. *et al.* (1981) Heavy chain variable region contribution to the NPb family of antibodies: somatic mutation evident in a gamma 2a variable region. *Cell* **24**, 625–37.

## VII. References (continued)

6. Senapathy, P., Shapiro, M.B. and Harris, N.L. (1990) Splice junctions, branch point sites, and exons: sequence statistics, identification, and applications to genome project. *Meth. Enzymol.* **183**, 252–78.
7. Gross, M.K., Kainz, M.S. and Merrill, G.F. (1987) Introns are inconsequential to efficient formation of cellular thymidine kinase mRNA in mouse L cells. *Mol. Cell. Biol.* **7**, 4576–81.
8. Buchman, A.R. and Berg, P. (1988) Comparison of intron-dependent and intron-independent gene expression. *Mol. Cell. Biol.* **8**, 4395–405.
9. Evans, M.J. and Scarpulla, R.C. (1989) Introns in the 3′-untranslated region can inhibit chimeric CAT and beta-galactosidase gene expression. *Gene* **84**, 135–42.
10. Huang, M.T.F. and Gorman, C.M. (1990) Intervening sequences increase efficiency of RNA 3′ processing and accumulation of cytoplasmic RNA. *Nucl. Acids Res.* **18**, 937–47.
11. Huang, M.T.F. and Gorman, C.M. (1990) The simian virus 40 small-t intron, present in many common expression vectors, leads to aberrant splicing. *Mol. Cell. Biol.* **10**, 1805–10.
12. Proudfoot, N.J. (1991) Poly(A) signals. *Cell* **64**, 671–4.
13. Bernstein, P. and Ross, J. (1989) Poly(A), poly(A) binding protein and the regulation of mRNA stability. *Trends Biochem. Sci.* **14**, 373–7.
14. Jackson, R.J. and Standart, N. (1990) Do the poly(A) tail and 3′ untranslated region control mRNA translation? *Cell* **62**, 15–24.
15. Carswell, S. and Alwine, J.C. (1989) Efficiency of utilization of the simian virus 40 late polyadenylation site: effects of upstream sequences. *Mol. Cell. Biol.* **9**, 4248–58.
16. Farr, A. and Roman, A. (1991) A pitfall of using a second plasmid to determine transfection efficiency. *Nucl. Acids Res.* **20**, 920.
17. Behr, J.P. *et al.* (1989) Efficient gene transfer into mammalian primary endocrine cells with lipopolyamine-coated DNA. *Proc. Natl. Acad. Sci. USA* **86**, 6982–6.
18. Loeffler, J.P. *et al.* (1990) Lipopolyamine-mediated transfection allows gene expression studies in primary neuronal cells. *J. Neurochem.* **54**, 1812–15.
19. Graham, F.L. and van der Eb, A.J. (1973) A new technique for the assay of infectivity of human adenovirus 5 DNA. *Virology* **52**, 456–67.
20. Wigler, M. *et al.* (1977) Transfer of purified herpes virus thymidine kinase gene to cultured mouse cells. *Cell* **11**, 223–32.
21. McCutchan, J.H. and Pagano, J.S. (1968) Enhancement of the infectivity of simian virus 40 deoxyribonucleic acid with diethylaminoethyl-dextran. *J. Natl. Cancer Inst.* **41**, 351.
22. Al-Moslih, M.I. and Dubes, G.R. (1973) The kinetics of DEAE-dextran-induced cell sensitization to transfection. *J. Gen. Virol.* **18**, 189–93.
23. Luthman, H. and Magnusson, G. (1983) High efficiency polyoma DNA transfection of chloroquine treated cells. *Nucl. Acids Res.* **11**, 1295–308.
24. Kawai, S. and Nishizawa, M. (1984) New procedure for DNA transfection with polycation and dimethyl sulfoxide. *Mol. Cell. Biol.* **4**, 1172–4.
25. Aubin, R.J., Weinfeld, M. and Paterson, M.C. (1988) Factors influencing efficiency and reproducibility of polybrene-assisted gene transfer. *Som. Cell Mol. Genet.* **14**, 155–67.
26. Andreason, G.L. and Evans, G.A. (1988) Introduction and expression of DNA molecules in eukaryotic cells by electroporation. *BioTechniques* **6**, 650–60.
27. Neumann, E. *et al.* (1982) Gene transfer into mouse lymphoma cells by electroporation in high electric fields. *EMBO J.* **1**, 841–5.

- (a) Certain applications of this product may require licenses from others.
- (b) Licensed under U.S. Pat. Nos. 5,292,658, 5,418,155 and other patents.
- (c) U.S. Pat. No. 5,283,179, Australian Pat. No. 649289 and other patents. Certain applications of this product may require licenses from others.
- (d) U.S. Pat. No. 5,744,320.
- (e) U.S. Pat. Nos. 5,492,817, 5,665,563, Australian Pat. No. 660329 and other patents.
- (f) The method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.
- (g) U.S. Pat. Nos. 4,966,964, 5,019,556 and 5,266,687, which claim vectors encoding a portion of human placental ribonuclease inhibitor, are exclusively licensed to Promega Corporation.
- (h) U.S. Pat. No. 5,552,302, Australian Pat. No. 646803 and other patents.
- (i) The cationic lipid component of the TransFast™ Transfection Reagent is covered by U.S. Pat. Nos. 5,824,812, 5,869,715 and pending foreign patents.
- (j) The cationic lipid component of the Tfx™ Reagents is covered by U.S. Pat. Nos. 5,527,928, 5,744,625 and 5,892,071, Australian Pat. No. 704189 and other pending foreign patents.
- (k) Transfectam is a registered trademark of BioSeptra, S.A., the holder of a license from CNRS-ULP Strasbourg under U.S. Pat. No. 5,171,678 to sell the Transfectam® product for research purposes only. The Transfectam® product was developed by J.P. Behr and J.P. Loeffler and is covered by the aforementioned patent.
- (l) U.S. Pat. No. 5,391,487.
- (m) The CMV promoter and its use are covered under U.S. Pat. Nos. 5,168,062 and 5,385,839 owned by the University of Iowa Research Foundation, Iowa City, Iowa, and licensed FOR RESEARCH USE ONLY. Commercial users must obtain a license to these patents directly from the University of Iowa Research Foundation.
- (n) U.S. Pat. No. 5,670,356.

© 1996–2001 Promega Corporation. All Rights Reserved.

Dual-Luciferase, ProFection and TNT are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office. Tfx and TransFast are trademarks of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the U.S. Department of Health and Human Services.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.



<b>Promega Corporation</b>	
2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Fax	608-277-2516
Internet	www.promega.com