



pUB6/V5-His A, B, and C

Catalog no. V250-01 and V250-20

Rev. date: 30 November 2009

Manual part no. 25-0251

MAN000089

User Manual

Table of Contents

Kit Contents and Storage.....	iv
Introduction	1
Product Overview	1
Methods	2
Cloning into pUB6/V5-His A, B, and C	2
Transfection and Analysis.....	6
Creating Stable Cell Lines	8
Appendix	11
Human UbC Promoter.....	11
pUB6/V5-His Vector.....	12
pUB6/V5-His/ <i>lacZ</i>	14
Blasticidin	15
Recipe.....	16
Accessory Products	17
Technical Support.....	18
Purchaser Notification.....	19
References.....	20

Kit Contents and Storage

Shipping and Storage

pUB6/V5-His A, B, and C vectors are shipped on wet ice. Upon receipt, store vectors at -20°C .

Kit Contents

The pUB6/V5-His A, B, and C vectors and the pUB6/V5-His /*lacZ* control plasmid are supplied with each product (Cat. no. V250-01, V250-20). The pUB6/V5-His with Blasticidin Kit (Cat. no. V250-01) also includes blasticidin antibiotic. The amount of each reagent provided is listed below.

Cat. no.	Vector	Composition	Amount
V250-01	pUB6/V5-His A, B, and C	40 μL of 0.5 $\mu\text{g}/\mu\text{L}$ vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 μg
	pUB6/V5-His / <i>lacZ</i>	40 μL of 0.5 $\mu\text{g}/\mu\text{L}$ vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 μg
	Blasticidin	powder	50 mg
V250-20	pUB6/V5-His A, B, and C	40 μL of 0.5 $\mu\text{g}/\mu\text{L}$ vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 μg
	pUB6/V5-His / <i>lacZ</i>	40 μL of 0.5 $\mu\text{g}/\mu\text{L}$ vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 μg

Intended Use

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

Introduction

Product Overview

Description of the System

pUB6/V5-His A, B, and C are 5.5 kb vectors designed for overproduction of recombinant proteins in mammalian cell lines. Features of the vectors allow purification and detection of expressed proteins (see pages 12–13 for more information). High-level stable and transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- Human ubiquitin C promoter (hUbC) for high-level expression across a broad range of species and cell types (Schorpp *et al.*, 1996; Wulff *et al.*, 1990) (see page 11 for more information).
- Three reading frames to facilitate in-frame cloning with a C-terminal peptide encoding the V5 epitope and a polyhistidine (6×His) metal-binding tag.
- Blasticidin resistance gene (*bsd*) for selection of stable cell lines (see page 15 for more information).
- Episomal replication in cell lines that are latently infected with SV40 or that express the SV40 large T antigen (*e.g.* COS-1, COS-7).

The control plasmid, pUB6/V5-His/*lacZ*, is included for use as a positive control for transfection, expression, and detection in the cell line of choice.

Experimental Outline

Use the following outline to clone and express your gene of interest in pUB6/V5-His.

1. Consult the multiple cloning sites described on pages 3–5 to determine which vector (A, B, or C) should be used to clone your gene in frame with the C-terminal V5 epitope and polyhistidine tag.
 2. Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on 50 to 100 µg/mL ampicillin or 50 µg/mL blasticidin.
 3. Analyze your transformants for the presence of insert by restriction digestion.
 4. Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in frame with the C-terminal peptide.
 5. Transfect your construct into the cell line of choice. Generate a stable cell line, if desired.
 6. Test for expression of your recombinant gene by western blot analysis or other functional assay. For antibodies to the V5 epitope or the polyhistidine, C-terminal tag, see the page 17.
 7. To purify your recombinant protein, you may use metal-chelating resin such as ProBond™. ProBond™ resin is available separately (see page 17 for ordering information).
-

Methods

Cloning into pUB6/V5-His A, B, and C

General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

E. coli Strain for Transformation

Many *E. coli* strains are suitable for the propagation of this vector, including TOP10F', JM109, and INV α F'. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*) and endonuclease A deficient (*endA*).

For your convenience, TOP10F' is available as chemically competent or electrocompetent cells from Invitrogen (see page 17).

Transformation Method

You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids.

Maintaining pUB6/V5-His

To propagate and maintain the pUB6/V5-His vectors, use a small amount of the supplied 0.5 $\mu\text{g}/\mu\text{L}$ stock solution in TE, pH 8.0 to transform a *recA*, *endA* *E. coli* strain like TOP10F', DH5 α , JM109, or equivalent. Select transformants on LB plates containing 50–100 $\mu\text{g}/\text{mL}$ ampicillin or 50 $\mu\text{g}/\text{mL}$ blasticidin. Be sure to prepare a glycerol stock of plasmid-containing *E. coli* strain for long-term storage (see page 5).

Cloning Considerations

Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak 1990). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position -3 and the G at position +4 (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined.

(G/A)NNATGG

If you wish to express your protein WITHOUT the C-terminal peptide, be sure to include a stop codon.

Continued on next page

Cloning into pUB6/V5-His A, B, and C, Continued

Multiple Cloning Site of pUB6/V5-His A

Below is the multiple cloning site for pUB6/V5-His A. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. **Note that there is a stop codon between the *Bam*H I site and the *Bst*X I site.** The multiple cloning site has been confirmed by sequencing and functional testing. The vector sequence may be downloaded from www.invitrogen.com or from **Technical Support** (see page 18). For more information on the hUbC promoter, see page 11.

```

1121      TTTAGGCACC TTTTGAAATG TAATCATTG GGTCAATATG TAATTTTCAG TGTTAGACTA GTAAATTGTC CGCTAAATTC
                                     UB forward priming site
1201      TGGCCGTTTT TGGCTTTTTT GTTAGAC GAA GCT TGG TAC CGA GCT CGG ATC CAC TAG TCC AGT GTG GTG
                                     Hind III  Asp718 I  Kpn I      Sac I BamH I                               BstX I*
                                     Glu Ala Trp Tyr Arg Ala Arg Ile His *** Ser Ser Val Val
1270      GAA TTC TGC AGA TAT CCA GCA CAG TGG CGG CCG CTC GAG TCT AGA GGG CCC TTC GAA GGT AAG CCT
EcoR I      Pst I  EcoR V      BstX I*  Not I      Xho I      Xba I      Apa I  BstB I
Glu Phe Cys Arg Tyr Pro Ala Gln Trp Arg Pro Leu Glu Ser Arg Gly Pro Phe Glu Gly Lys Pro
1336      ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC CAT CAC CAT TGA GT
V5 epitope      Age I      Polyhistidine region
Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His His His His ***
1401      TTAAACCCGC TGATCAGCCT CGACTGTGCC TTCTAGTTGC CAGCCATCTG TTGTTTGCCC CTCCCCCGTG CCTTCCTTGA
Pme I      BGH reverse priming site
1481      CCCTGGAAGG TGCCACTCCC ACTGTCCTTT CTAATAAAA TGAGGAAATT GCATCGCATT GTCTGAGTAG GTGTCATTCT
                                     BGH polyadenylation signal

```

*Note that there are two *Bst*X I sites in the polylinker.

Continued on next page

Cloning into pUB6/V5-His A, B, and C, Continued

Multiple Cloning Site of pUB6/V5-His B

Below is the multiple cloning site for pUB6/V5-His B. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. The multiple cloning site has been confirmed by sequencing and functional testing. The vector sequence may be downloaded from www.invitrogen.com or from **Technical Support** (see page 18). For more information on the hUbC promoter, see page 11.

```

1121  TTTAGGCACC TTTTCAAATG TAATCATTG GGTCAATATG TAATTTTCAG TGTTAGACTA GTAAATTGTC CGCTAAATTC
                                     UB forward priming site
1201  TGGCCGTTTT TGGCTTTTTT GTTAGACG AAG CTT GGT ACC GAG CTC GGA TCC ACT AGT CCA GTG TGG TGG
                                     Hind III  Asp718 I  Kpn I      Sac I BamH I                               BstX I* EcoR I
                                     Lys Leu Gly Thr Glu Leu Gly Ser Thr Ser Pro Val Trp Trp
1271  AAT TCT GCA GAT ATC CAG CAC AGT GGC GGC CGC TCG AGT CTA GAG GGC CCG CCG TTC GAA GGT AAG
      Pst I  EcoRV      BstX I*  Not I      Xho I      Xba I      Apa I      BstB I
      Asn Ser Ala Asp Ile Gln His Ser Gly Gly Arg Ser Ser Leu Glu Gly Pro Arg Phe Glu Gly Lys
1337  CCT ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC CAT CAC CAT TGA
      V5 epitope                                     Age I      Polyhistidine region
      Pro Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His His His His ***
1403  GTTTAAAC CCGCTGATCA GCCTCGACTG TGCCTTCTAG TTGCCAGCCA TCTGTTGTTT GCCCCTCCCC CGTGCCTTCC
      Pme I      BGH reverse priming site
1481  TTGACCCTGG AAGGTGCCAC TCCCACTGTC CTTTCCTAAT AAAATGAGGA AATTGCATCG CATTGTCTGA GTAGGTGTC
                                     BGH polyadenylation signal

```

Note that there are two *BstX I* sites in the polylinker.

Continued on next page

Cloning into pUB6/V5-His A, B, and C, Continued

Multiple Cloning Site of pUB6/V5-His C

Below is the multiple cloning site for pUB6/V5-His C. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. The multiple cloning site has been confirmed by sequencing and functional testing. The vector sequence may be downloaded from www.invitrogen.com or from **Technical Support** (see page 18).

```

1121      TTTAGGCACC TTTTGAATG TAATCATTG GGTC AATATG TAATTTTCAG TGTTAGACTA GTAAATTGTC CGCTAAATTC
                                     UB forward priming site
1201      TGGCCGTTTT TGGCTTTTTT GTTAGACGA AGC TTG GTA CCG AGC TCG GAT CCA CTA GTC CAG TGT GGT GGA
                                     Hind III  Asp718 I  Kpn I      Sac I  BamH I                                     BstX I*  EcoR I
                                     Ser Leu Val Pro Ser Ser Asp Pro Leu Val Gln Cys Gly Gly
1272      ATT CTG CAG ATA TCC AGC ACA GTG GCG GCC GCT CGA GGT CAC CCA TTC GAA GGT AAG CCT ATC CCT
                                     Pst I  EcoR V      BstX I*  Not I      Xho I      BstE II      BstB I      V5 epitope
                                     Ile Leu Gln Ile Ser Ser Thr Val Ala Ala Ala Arg Gly His Pro Phe Glu Gly Lys Pro Ile Pro
1338      AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC CAT CAC CAT TGA GTTTAA
                                     Age I      Polyhistidine region                                     Pme I
                                     Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His His His His ***
1401      ACCCGCTGAT CAGCCTCGAC TGTGCCTTCT AGTTGCCAGC CATCTGTTGT TTGCCCTCC CCCGTGCCTT CCTTGACCTT
                                     BGH Reverse priming site
1481      GGAAGGTGCC ACTCCCACTG TCCTTTCCTA ATAAAATGAG GAAATTGCAT CGCATTGTCT GAGTAGGTGT CATTCTATTC
                                     BGH polyadenylation signal

```

*Note that there are two *BstX I* sites in the polylinker.

Transforming Ligation Mixtures

Transform your ligation mixtures into a competent *recA, endA E. coli* strain (e.g., TOP10F⁺, DH5 α) and select on LB plates containing 50–100 $\mu\text{g}/\text{mL}$ ampicillin or 50 $\mu\text{g}/\text{mL}$ blasticidin. Select 10–20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the UB Forward and BGH Reverse primer sequences to confirm that your gene is fused in frame with the V5 epitope and the C-terminal polyhistidine tag. Refer to the diagram on pages 3–5 for the sequence and location of the primer binding sites.

Preparing a Glycerol Stock

Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. It is also a good idea to keep a DNA stock of your plasmid at -20°C .

1. Streak the original colony out for single colonies on an LB plate containing 50 $\mu\text{g}/\text{mL}$ ampicillin (or 50 $\mu\text{g}/\text{mL}$ blasticidin). Incubate the plate at 37°C overnight.
2. Isolate a single colony and inoculate into 1–2 mL of LB containing 50 $\mu\text{g}/\text{mL}$ ampicillin (or 50 $\mu\text{g}/\text{mL}$ blasticidin).
3. Grow the culture to mid-log phase ($\text{OD}_{600} = 0.5\text{--}0.7$).
4. Mix 0.85 mL of culture with 0.15 mL of sterile glycerol and transfer to a cryovial.
5. Store at -80°C .

Transfection and Analysis

Introduction

Once you have confirmed that your construct is in the correct orientation and fused to the C-terminal peptide (if desired), then you are ready to transfect your cell line of choice. We recommend that you include the positive control vector and a mock transfection to evaluate your results.

Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink™ HiPure Miniprep Kit or the PureLink™ HiPure Midiprep Kit (see page 17 for ordering information).

Methods of Transfection

For established cell lines (e.g., HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. It is recommended that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in Current Protocols in Molecular Biology.

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler et al., 1977), lipid-mediated (Felgner et al., 1987; Felgner et al., 1989), and electroporation (Chu et al., 1987; Shigekawa and Dower, 1988). Invitrogen offers the Lipofectamine™ 2000 Reagent for mammalian transfection (see page 17 for ordering).

Positive Control

pUB6/V5-His/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 14). It may be used to optimize transfection conditions for your cell line. The gene encoding β -galactosidase is expressed in mammalian cells under the control of the human ubiquitin C (hUbC) promoter. A successful transfection will result in β -galactosidase expression that can be easily assayed (see below).

Assay for β -galactosidase Activity

You may assay for β -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the β -Gal Assay Kit and the β -Gal Staining Kit for fast and easy detection of β -galactosidase expression (see page 17 for ordering information).

Continued on next page

Transfection and Analysis, Continued

Detecting Fusion Proteins

Several antibodies are available from Invitrogen that can be used to detect expression of your fusion protein from pUB6/V5-His (see page 17 for ordering information).

To detect the fusion protein by western blot, you will need to prepare a cell lysate from transfected cells. We recommend that you perform a time course to optimize expression of the fusion protein (*e.g.* 24, 48, 72 hours, etc. after transfection). To lyse cells:

1. Wash cell monolayers ($\sim 10^6$ cells) once with phosphate-buffered saline (PBS).
 2. Scrape cells into 1 mL PBS and pellet the cells at $1,500 \times g$ for 5 minutes.
 3. Resuspend in 50 μ L Cell Lysis Buffer (see recipe below). Other lysis buffers are suitable.
 4. Incubate cell suspension at 37°C for 10 minutes to lyse the cells.
 5. Vortex the cell lysate and centrifuge at $10,000 \times g$ for 10 minutes to pellet nuclei. Transfer the supernatant to a fresh tube. Assay the supernatant for protein concentration. **Note:** Do not use protein assays utilizing Coomassie[®] Blue or other dyes. NP-40 interferes with the binding of the dye with the protein.
 6. Add SDS-PAGE sample buffer to a final concentration of 1X and boil the sample for 5 minutes.
 7. Load 20 μ g of lysate onto an SDS-PAGE gel and electrophorese. Use the appropriate percentage of acrylamide to resolve your fusion protein.
-

Creating Stable Cell Lines

Introduction

The pUB6/V5-His vectors contain the blasticidin resistance gene (*bsd*) to allow for selection of stable cell lines using blasticidin (Kimura *et al.*, 1994). We recommend that you test the sensitivity of your mammalian host cell to blasticidin as natural resistance varies among cell lines. General information and guidelines are provided below for your convenience.

Possible Sites for Linearization

To obtain stable transfectants, you may choose to linearize your vector before transfection. While linearizing your vector may not improve the efficiency of transfection, it increases the chances that the vector does not integrate in a way that disrupts the gene of interest. The table below lists some unique sites that may be used to linearize your construct prior to transformation. **Other restriction sites are possible.** Note that for the enzymes listed below, the cleavage site is indicated for versions A, B, and C of pUB6/V5-His. Be sure that your insert does not contain the restriction enzyme site you wish to use to linearize your vector.

Enzyme	Location	Supplier
<i>Bgl</i> II	Upstream of hUbc promoter	Many
<i>Bst</i> 1107 I	End of SV40 polyA	AGS*, Fermentas, Takara
<i>Sap</i> I	Backbone	New England Biolabs
<i>Bsp</i> LU 11I	Backbone	Boehringer-Mannheim
<i>Alw</i> N I	pMB1 origin	Amersham, New England Biolabs, Life Technologies
<i>Eam</i> 1105 I	Ampicillin gene	AGS*, Fermentas, Takara
<i>Bgl</i> I	Ampicillin gene	Many
<i>Fsp</i> I	Ampicillin gene	Many
<i>Sca</i> I	Ampicillin gene	Many
<i>Ssp</i> I	Backbone	Many

Continued on next page

Creating Stable Cell Lines, Continued

Selection in Mammalian Cell Lines

To generate a stable cell line expressing your protein, you need to determine the minimum concentration of blasticidin required to kill your untransfected host cell line. Typically, concentrations between 2 and 10 $\mu\text{g}/\text{mL}$ blasticidin are sufficient to kill the untransfected host cell line. Test a range of concentrations (see below) to ensure that you determine the minimum concentration necessary for your cell line. See page 15 for details on handling and preparing Blasticidin solution.

1. Seed cells at 20–25% confluency for each time point (~6 time points) and allow the cells to adhere overnight.
 2. The next day, substitute culture medium with medium containing varying concentrations of blasticidin (*e.g.*, 0, 1, 3, 5, 7.5, and 10 $\mu\text{g}/\text{mL}$ blasticidin).
 3. Replenish the selective medium every 3–4 days. Cells sensitive to blasticidin will round up and detach from the plate. Dead cells will accumulate in the medium.
 4. Count the number of viable cells at regular intervals to determine the appropriate concentration of blasticidin that prevents growth.
-

Selecting Stable Integrants

Once the appropriate concentration of blasticidin is determined, you can generate a stable cell line with your construct.

1. Transfect your cells using the optimal protocol for your cell line. Include a sample of untransfected cells as a negative control.
 2. 48 hours after transfection, split the cells into fresh medium containing blasticidin at the appropriate concentration for your cell line. Split the cells such that they are no more than 25% confluent.
 3. Replenish selective medium every 3–4 days until blasticidin-resistant colonies are detected.
 4. Pick and expand colonies.
-

Continued on next page

Creating Stable Cell Lines, Continued

Preparing Cells for Lysis

Use the procedure below to prepare cells for lysis prior to purification of your protein on ProBond™. You will need 5×10^6 to 1×10^7 cells for purification of your protein on a 2 mL ProBond™ column (see ProBond™ Protein Purification manual).

1. Seed cells (from a stable cell line) in five T-75 flasks or 2 to 3 T-175 flasks.
 2. Grow the cells in selective medium until they are 80–90% confluent.
 3. Harvest the cells by treating with trypsin-EDTA for 2 to 5 minutes or by scraping the cells in PBS.
 4. Inactivate the trypsin by diluting with fresh medium (if necessary) and transfer the cells to a sterile microcentrifuge tube.
 5. Centrifuge the cells at $1,500 \times g$ for 5 minutes. Resuspend the cell pellet in PBS.
 6. Centrifuge the cells at $1,500 \times g$ for 5 minutes. You may lyse the cells immediately or freeze in liquid nitrogen and store at -80°C until needed.
-

Lysis of Cells

If you are using ProBond™ resin, refer to the ProBond™ Protein Purification manual for details about sample preparation for chromatography.


If you are using other metal-chelating resin, refer to the manufacturer's instruction for recommendations on sample preparation.

Appendix

Human UbC Promoter

Description


The human UbC promoter allows high-level expression of recombinant protein in most mammalian cell lines (Wulff *et al.*, 1990) and in virtually all tissues tested in transgenic mice (Schorpp *et al.*, 1996). The diagram below shows the features of the UbC promoter used in pUB6/V5-His (Nenoi *et al.*, 1996). Features are marked as per Nenoi, *et al.*, 1996.



pUB6/V5-His Vector

Map of pUB6/V5-His

The figure below summarizes the features of the pUB6/V5-His vectors. The sequences for pUB6/V5-His A, B, and C are available for downloading from www.invitrogen.com or from **Technical Support** (see page 18).



Comments for pUB6/V5-His A: 5463 nucleotides

- UbC promoter: bases 18-1227
- UB forward priming site: bases 1167-1188
- Multiple cloning site: bases 1229-1326
- V5 epitope: bases 1327-1368
- Polyhistidine tag: bases 1378-1395
- BGH reverse priming site: bases 1418-1435
- BGH polyadenylation signal: bases 1421-1648
- f1 origin: bases 1694-2122
- SV40 promoter and origin: bases 2149-2458
- EM-7 promoter: bases 2506-2561
- Blasticidin resistance gene: bases 2580-2978
- SV40 polyadenylation signal: bases 3136-3266
- pUC origin: bases 3649-4322
- Ampicillin resistance gene: bases 4467-5327

* Frame-dependent variations.
BstE II is only found in version C.
In addition, there are no *Xba I*
or *Apa I* sites in version C.

Continued on next page

pUB6/V5-His Vector, Continued

Features of pUB6/V5-His

pUB6/V5-His A (5463 bp), pUB6/V5-His B (5467 bp), and pUB6/V5-His C (5459 bp) contain the following elements. All features have been functionally tested.

Feature	Benefit
Human ubiquitin C (hUbC) promoter	Allows overexpression of your recombinant protein in a broad range of mammalian cell types (Hershko and Ciechanover, 1982; Wulff <i>et al.</i> , 1990; Schorpp <i>et al.</i> , 1996).
Multiple cloning site in three reading frames	Allows insertion of your gene and facilitates cloning in frame with the V5 epitope and the C-terminal polyhistidine tag.
V5 epitope (Gly-Lys-Pro-Ile-Pro-Asn-Pro-Leu-Leu-Gly-Leu-Asp-Ser-Thr)	Allows detection of your recombinant protein with the Anti-V5 Antibodies (Southern <i>et al.</i> , 1991).
C-terminal polyhistidine tag	Allows purification of your recombinant protein on metal-chelating resin such as ProBond™. In addition, the C-terminal polyhistidine tag is the epitope for the Anti-His(C-term) Antibodies (Lindner <i>et al.</i> , 1997).
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992).
f1 origin	Allows rescue of single-stranded DNA.
SV40 early promoter and origin	Allows efficient, high-level expression of the blasticidin resistance gene and episomal replication in cells expressing the SV40 large T antigen.
EM-7 promoter	For expression of the blasticidin resistance gene in <i>E. coli</i> .
Blasticidin resistance gene (<i>bsd</i>)	Selection of stable transfectants in mammalian cells (Kimura <i>et al.</i> , 1994).
SV40 polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA.
pUC origin	High-copy number replication and growth in <i>E. coli</i> .
Ampicillin resistance gene (β-lactamase)	Selection of vector in <i>E. coli</i> .


pUB6/V5-His/lacZ

Description

pUB6/V5-His/lacZ is a 8510 bp control vector containing the gene for β -galactosidase. This vector was constructed by ligating a 3190 bp *Hind* III-*Age* I fragment containing the *lacZ* gene and the V5 epitope to a 5320 bp *Hind* III-*Age* I fragment containing the hUbC promoter, polyhistidine tag and blasticidin resistance gene from pUB6/V5-His A.

Map of Control Vector

The figure below summarizes the features of the pUB6/V5-His/lacZ vector. The nucleotide sequence for pUB6/V5-His/lacZ is available for downloading from www.invitrogen.com or by contacting **Technical Support**. (see page 18).



Comments for pUB6/V5-His/lacZ: 8510 nucleotides

- UbC promoter: bases 18-1227
- UB forward priming site: bases 1167-1188
- LacZ* ORF: bases 1290-4346
- V5 epitope: bases 4374-4415
- Polyhistidine tag: bases 4425-4442
- BGH reverse priming site: bases 4465-4482
- BGH polyadenylation signal: bases 4468-4695
- r1 origin: bases 4741-5169
- SV40 promoter and origin: bases 5196-5505
- EM-7 promoter: bases 5553-5608
- Blasticidin resistance gene: bases 5627-6025
- SV40 polyadenylation signal: bases 6183-6313
- pUC origin: bases 6696-7369
- Ampicillin resistance gene: bases 7514-8374

Blasticidin

Description

Blasticidin S HCl is a nucleoside antibiotic isolated from *Streptomyces griseochromogenes* which inhibits protein synthesis in both prokaryotic and eukaryotic cells (Takeuchi *et al.*, 1958; Yamaguchi *et al.*, 1965). Resistance is conferred by expression of either one of two blasticidin S deaminase genes: *BSD* from *Aspergillus terreus* (Kimura *et al.*, 1994) or *bsr* from *Bacillus cereus* (Izumi *et al.*, 1991). These deaminases convert blasticidin S to a non-toxic deaminohydroxy derivative (Izumi *et al.*, 1991).

Handling Blasticidin

Always wear gloves, mask, goggles, and protective clothing (e.g., a laboratory coat) when handling blasticidin. Weigh out blasticidin and prepare solutions in a hood.


Preparing and Storing Stock Solutions

Blasticidin is soluble in water. Water is generally used to prepare stock solutions of 5 to 10 mg/mL.

- Dissolve blasticidin in sterile water and filter-sterilize the solution.
 - Aliquot in small volumes suitable for one time use (see last point below) and freeze at -20°C for long-term storage or store at 4°C for short term storage.
 - Aqueous stock solutions are stable for 1–2 weeks at 4°C and 6–8 weeks at -20°C .
 - pH of the aqueous solution should not exceed 7 to prevent inactivation of blasticidin.
 - Do not subject stock solutions to freeze/thaw cycles (do not store in a frost-free freezer).
 - Upon thawing, use what you need and discard the unused portion.
-

Molecular Weight, Formula, and Structure

The formula for blasticidin is $\text{C}_{17}\text{H}_{26}\text{N}_8\text{O}_5\text{-HCl}$, and the molecular weight is 458.9. The diagram below shows the structure of blasticidin.



Recipe

Cell Lysis Buffer

50 mM Tris
150 mM NaCl
1% Nonidet P-40
pH 7.8

1. This solution can be prepared from the following common stock solutions.
For 100 mL, combine:
1 M Tris base 5 mL
5 M NaCl 3 mL
Nonidet P-40 1 mL
2. Bring the volume to 90 mL with deionized water and adjust the pH to 7.8 with HCl.
3. Bring the volume up to 100 mL. Store at room temperature.

Note: Protease inhibitors may be added at the following concentrations:

1 mM PMSF
1 µg/mL Pepstatin
1 µg/mL Leupeptin

Accessory Products

Introduction

The following products may be used with the pUB6/V5-His vectors. For details, visit www.invitrogen.com or contact **Technical Support** (page 18).

Item	Amount	Catalog no.
ProBond™ Purification System	6 × 2 mL precharged, prepacked ProBond™ resin columns and buffers for native and denaturing purification	K850-01
ProBond™ Resin	50 mL	R801-01
	150 mL	R801-15
PureLink™ HiPure Plasmid Midiprep Kit	25 preps	K2100-04
Lipofectamine™ 2000 Reagent	0.75 mL	11668-027
Electrocomp™ TOP10F'	5 × 80 µL	C665-55
One Shot® TOP10F' Chemically Competent <i>E. coli</i>	20×50 µL	C3030-03
β-Gal Assay Kit	80 mL	K1455-01
β-Gal Staining Kit	1 kit	K1465-01

Antibodies

If you do not have an antibody specific to your protein, Invitrogen offers the Anti-V5, or Anti-His(C-term) antibodies to detect your recombinant fusion protein. Horseradish peroxidase (HRP)- and alkaline phosphatase (AP)-conjugated antibodies are available for convenient one-step detection.

Antibody	Epitope	Catalog no.
Anti-V5	Detects a 14 amino acid epitope derived from the P and V proteins of the paramyxovirus, SV5 (Southern et al., 1991): GKPIP NPLLGLDST	R960-25
Anti-V5-HRP		R961-25
Anti-V5-AP		R962-25
Anti-His(C-term)	Detects the C-terminal polyhistidine tag (requires the free carboxyl group for detection) (Lindner et al., 1997): HHHHHH-COOH	R930-25
Anti-His(C-term)-HRP		R931-25
Anti-His(C-term)-AP		R932-25

Primers

For your convenience, Invitrogen offers a custom primer synthesis service. Visit www.invitrogen.com for more details.

Technical Support

Web Resources



Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
 - Complete technical support contact information
 - Access to the Invitrogen Online Catalog
 - Additional product information and special offers
-

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website (www.invitrogen.com).

Corporate Headquarters:

5791 Van Allen Way
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail: tech_support@invitrogen.com

Japanese Headquarters:

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Tel: 81 3 5730 6509
Fax: 81 3 5730 6519
E-mail: jpinfo@invitrogen.com

European Headquarters:

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

SDS

Safety Data Sheets (MSDSs) are available on our website at www.invitrogen.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.invitrogen.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives. All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

**Limited Use Label
License No. 22:
Vectors and
Clones Encoding
Histidine Hexamer**

This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

**Limited Use Label
License No. 51:
Blasticidin and the
Blasticidin
Selection Marker**

Blasticidin and the blasticidin resistance gene (bsd) are the subject of U.S. Patent No.5,527,701 sold under patent license for research purposes only. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. email: outlicensing@invitrogen.com

References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nuc. Acids Res.* *15*, 1311-1326.
- Felgner, P. L., Gadek, T. R., Holm, M., Roman, R., Chan, H. W., Wenz, M., Northrop, J. P., Ringold, G. M., and Danielsen, M. (1987). Lipofectin: A Highly Efficient, Lipid-mediated DNA-transfection Procedure. *Proc. Natl. Acad. Sci. USA* *84*, 7413-7417.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Hershko, A., and Ciechanover, A. (1982). Mechanisms of Intracellular Protein Breakdown. *Ann. Rev. Biochem.* *51*, 335-364.
- Izumi, M., Miyazawa, H., Kamakura, T., Yamaguchi, I., Endo, T., and Hanaoka, F. (1991). Blastocidin S-Resistance Gene (*bsr*): A Novel Selectable Marker for Mammalian Cells. *Exp. Cell Res.* *197*, 229-233.
- Kimura, M., Takatsuki, A., and Yamaguchi, I. (1994). Blastocidin S Deaminase Gene from *Aspergillus terreus* (*BSD*): A New Drug Resistance Gene for Transfection of Mammalian Cells. *Biochim. Biophys. Acta* *1219*, 653-659.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nuc. Acids Res.* *15*, 8125-8148.
- Kozak, M. (1991). An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biol.* *115*, 887-903.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Lindner, P., Bauer, K., Krebber, A., Nieba, L., Kremmer, E., Krebber, C., Honegger, A., Klinger, B., Mocikat, R., and Pluckthun, A. (1997). Specific Detection of His-tagged Proteins With Recombinant Anti-His Tag scFv-Phosphatase or scFv-Phage Fusions. *BioTechniques* *22*, 140-149.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nenoi, M., Mita, K., Ichimura, S., Cartwright, I. L., Takahashi, E., Yamaguchi, M., and Tsuji, H. (1996). Heterogeneous Structure of the Polyubiquitin Gene UbC of HeLa S3 Cells. *Gene* *175*, 179-185.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Schorpp, M., Jäger, R., Schellander, K., Schenkel, J., Wagner, E. F., Weiher, H., and Angel, P. (1996). The Human Ubiquitin C Promoter Directs High Ubiquitous Expression of Transgenes in Mice. *Nuc. Acids Res.* *24*, 1787-1788.

Continued on next page

References, Continued

- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* 6, 742-751.
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991). Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. *J. Gen. Virol.* 72, 1551-1557.
- Takeuchi, S., Hirayama, K., Ueda, K., Sakai, H., and Yonehara, H. (1958). Blasticidin S, A New Antibiotic. *The Journal of Antibiotics, Series A* 11, 1-5.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* 11, 223-232.
- Wulff, B. S., O'Hare, M. M., Boel, E., Theill, L. E., and Schwartz, T. W. (1990). Partial Processing of the Neuropeptide Y Precursor in Transfected CHO Cells. *FEBS Lett.* 261, 101-105.
- Yamaguchi, H., Yamamoto, C., and Tanaka, N. (1965). Inhibition of Protein Synthesis by Blasticidin S. I. Studies with Cell-free Systems from Bacterial and Mammalian Cells. *J. Biochem. (Tokyo)* 57, 667-677.

©2009 Life Technologies Corporation. All rights reserved.

Notes



Corporate Headquarters

5791 Van Allen Way

Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information, visit our web site at www.invitrogen.com