

Use

pLEGFP-C1 or derivatives expressing EGFP fusions to your gene of interest can be transfected into any of BD Biosciences Clontech packaging cell lines (see www.bdbiosciences.com for a complete list). pLEGFP-C1 transcripts produced in the packaging cell contain the Ψ^+ (psi) RNA packaging signal, the neomycin gene, transcription and processing elements, and the gene of interest. pLEGFP-C1 does not contain the structural *gag*, *pol*, and *env* genes necessary for retroviral particle formation and replication; however, these genes are stably integrated in the packaging cell genome. Once in the cell, RNA from the vector is packaged into high-titer, infectious, replication-incompetent retroviral particles (9–12). That is, these retroviral particles can infect target cells and transmit the EGFP or EGFP-fusion gene, but cannot replicate within these cells since the cells lack the viral structural genes. The separate introduction and integration of the structural genes into the packaging cell line minimizes the chances of producing replication-competent virus due to recombination events during cell proliferation.

Location of features

- 5' MoMuLV LTR: 145–733
- Ψ^+ (extended packaging signal): 803–1612
- Neomycin resistance gene (Neo^r): 1656–2450
- Immediate early CMV promoter (P_{CMV}): 2468–3057
- Enhanced green fluorescent protein (EGFP) gene: 3074–3790
 - Kozak consensus translation initiation site: 2931–2941
 - Start codon (ATG): 3074–3076;
 - Last EGFP codon: 3788–3790
 - Stop codon of the protein expressed by the native vector: 3869–3871
 - (Note that the version of EGFP encoded by this vector differs slightly from the original EGFP due to the extra amino-acids at the C-terminus)
 - Stop codons downstream of the MCS: 3865–3867; 3869–3871; 3873–3875;
 - (Note that stop codons may be introduced by inserts)
 - Insertion of Val at position 2: 3077–3079
 - GFPmut1 chromophore mutations (Phe-64 to Leu; Ser-65 to Thr): 3355–3360
 - His-231 to Leu mutation (A→T): 3768
- Multiple Cloning Site: 3799–3868
- 3' MoMuLV LTR: 3938–4531
- pBR322 plasmid replication region
 - Site of replication initiation: 5068
- Ampicillin resistance gene (β -lactamase): 6687–5827

Sequencing primer locations

- EGFP-C Sequencing Primer (#6478-1 [3727–3748]):
5'-CATGGTCCTGCTGGAGTTCGTG-3'
- 3' primer pLNCX Seq/PCR Primer (#K1060-F; [3960–3935]):
5'-ACCTACAGGTGGGGTCTTTCATTCCC-3'
(The 5' primer in this set anneals upstream of the EGFP sequence and is not useful for analyzing derivatives of pLEGFP-C1.)

Propagation in *E. coli*

- Suitable host strains: DH5 α , HB101, and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 μ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: pBR322
- Copy number: low

NOTE: The viral supernatants produced by this retroviral vector could, depending on your cloned insert, contain potentially hazardous recombinant virus. Due caution must be exercised in the production and handling of recombinant retrovirus. Appropriate NIH, regional, and institutional guidelines apply.

References:

1. Cormack, B., *et al.* (1996) *Gene* **173**:33–38.
2. Haas, J., *et al.* (1996) *Curr. Biol.* **6**:315–324.
3. Kozak, M. (1987) *Nucleic Acids Res.* **15**:8125–8148.
4. *Clontechniques* (April 1996) **XI** (2):2–3.
5. Coffin, J. M. & Varmus, H. E., Eds. (1996) *Retroviruses* (Cold Spring Harbor Laboratory Press, NY).
6. Ausubel, F. M., *et al.* (1994) *Current Protocols in Molecular Biology* (Greene Publishing Associates, Inc. & John Wiley & Sons, Inc.).
7. Miller, A. D. & Rosman, G. J. (1989) *BioTechniques* **7**:980–990.
8. *Clontechniques* (April 1997) **XII** (2):8–9.
9. Mann, R., *et al.* (1983) *Cell* **33**:153–159.
10. Miller, A. D. & Buttimore, C. (1986) *Mol. Cell. Biol.* **6**:2895–2902.
11. Morgenstern, J. P. & Land, H. (1990) *Nucleic Acids Res.* **18**:3587–3590.
12. Miller, A. D. & Chen, F. (1996) *J. Virol.* **70**:5564–5571.

Notice to Purchaser

Use of BD Biosciences Clontech's Living Colors™ products containing DNA sequences coding for mutant *Aequorea victoria* green fluorescent protein (GFP) variants or proteins thereof requires a license from Amersham Biosciences under U.S. Patent Nos. 5,625,048; 5,777,079; 6,054,321 and other pending U.S. and foreign patent applications. In addition, certain BD Biosciences Clontech products are made under U.S. Patent No. 5,804,387 licensed from Stanford University.

Not-For-Profit research institutes or entities are granted an automatic license with the purchase of this product for use in non-commercial internal research purposes, the terms of which are disclosed in detail in the license that accompanies the shipment of this product. Such license specifically excludes the right to sell or otherwise transfer this product or its components to third parties.

For-Profit research institutes or entities must obtain a license from Amersham Biosciences. E-mail: gfp@amershambiosciences.com

Please contact BD Biosciences Clontech directly for any other assistance, including purchasing and technical support. All companies and institutions purchasing Living Colors™ products will be included in a quarterly report to Aurora Biosciences, as required by the BD Biosciences Clontech/Aurora Biosciences license agreement.

The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by BD Biosciences Clontech. This vector has not been completely sequenced.

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. BD Biosciences Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of BD Biosciences Clontech.

© 2002, Becton, Dickinson and Company